APPLICATION FOR UNITED STATES LETTERS PATENT

(Case no. 03-389-A)

5

TITLE:	METHOD AND SYSTEM FOR COMPREHENSIVE EVALUATION
1 1 1 1 1 1	

OF ORTHODONTIC TREATMENT USING UNIFIED

WORKSTATION

10

INVENTORS: Rohit Sachdeva, a citizen of the United States and a resident of Plano,

1

Texas

Peer Sporbert, a citizen of Germany and a resident of Berlin, Germany

15

Stephan Maetzel, a citizen of Germany and a resident of Berlin,

Germany

Hans Imgrund, a citizen of Germany and a resident of Berlin, Germany

20

Phillip Getto, a citizen of the United States and a resident of Plano,

Texas

Doke Evan Roberts, a citizen of the United States and a resident of

Lucas, Texas

30

25

Assignee: OraMetrix, Inc.

BACKGROUND OF THE INVENTION

Related applications

This is a continuation-in-part of application serial no. 10/428,461 filed May 2, 2003,

which is a continuation-in-part of application serial no. 09/834,412, filed April 13, 2001.

5 The entire contents of both related applications are incorporated by reference herein.

A. Field of the Invention

10

15

20

25

This invention relates to the field of computerized techniques for orthodontic

treatment planning for human patients. More particularly, the invention is directed to an

interactive workstation and associated computerized techniques for facilitating integration of

various tasks performed in planning treatment for orthodontic patients, and more particularly

evaluation of a proposed treatment plan.

B. Description of Related Art

The traditional process of diagnosis and treatment planning for a patient with

orthodontic problems or disease typically consists of the practitioner obtaining clinical

history, medical history, dental history, and orthodontic history of the patient supplemented

by 2D photographs, 2D radiographic images, CT scans, 2D and 3D scanned images,

ultrasonic scanned images, and in general non-invasive and sometimes invasive images, plus

video, audio, and a variety of communication records. Additionally, physical models, such as

made from plaster of paris, of the patient's teeth are created from the impressions taken of the

patient's upper and lower jaws. Such models are manually converted into teeth drawings by

projecting teeth on drawing paper. Thus, there is a large volume of images and data involved

in the diagnosis and treatment planning process. Furthermore, the information may require

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

conversion from one form to another and selective reduction before it could become useful.

There are some computerized tools available to aid the practitioner in these data conversion

and reduction steps, for example to convert cephalometric x-rays (i.e., 2 dimensional x-ray

photographs showing a lateral view of the head and jaws, including teeth) into points of

interest with respect to soft tissue, hard tissue, etc., but they are limited in their functionalities

and scope. Even then, there is a fairly substantial amount of manual work involved in these

steps.

10

15

20

Additionally, a number of measurements, e.g., available space between teeth, are also

often done manually. Generally, these steps are time consuming and prone to inherent

inaccuracies. Furthermore, the practitioner has to contend with the biological

interdependencies within the patient, which introduces constraints eliminating certain

treatment options that would otherwise be acceptable, between the soft tissue, the hard tissue,

and the teeth. There is lack of an integrated platform which a practitioner could utilize to

filter-out non-practicable treatment options.

Consequently, the practitioner is left to mental visualization, chance process to select

the treatment course that would supposedly work. Furthermore, the diagnosis process is

some-what ad-hoc and the effectiveness of the treatment depends heavily upon the

practitioner's level of experience. Often, due to the complexities of the detailed steps and the

time consuming nature of them, some practitioners take a short-cut, relying predominantly on

their intuition to select a treatment plan. For example, the diagnosis and treatment planning

is often done by the practitioner on a sheet of acetate over the X-rays. All of these factors

frequently contribute towards trial and error, hit-and-miss, lengthy and inefficient treatment

plans that require numerous mid-course adjustments. While at the beginning of treatment

things generally run well as all teeth start to move at least into the right direction, at the end

3

of treatment a lot of time is lost by adaptations and corrections required due to the fact that

the end result has not been properly planned at any point of time. By and large, this approach

lacks reliability, reproducibility and precision. More over, there is no comprehensive way

available to a practitioner to stage and simulate the treatment process in advance of the actual

implementation to avoid the often hidden pitfalls. And the patient has no choice and does not

know that treatment time could be significantly reduced if proper planning was done.

In recent years, computer-based approaches have been proposed for aiding

orthodontists in their practice. However, these approaches are limited to diagnosis and

treatment planning of craniofacial structures, including the straightening of teeth. See

Andreiko, U.S. Patent 6,015,289; Snow, U.S. Patent No. 6,068,482; Kopelmann et al., U.S.

Patent 6,099,314; Doyle, et al., U.S. Patent 5,879,158; Wu et al., U.S. Patent 5,338,198, and

Chisti et al., U.S. Patents 5,975,893 and 6,227,850, the contents of each of which is

incorporated by reference herein. Also see imaging and diagnostic software and other related

products marketed by Dolphin Imaging, 6641 Independence Avenue, Canoga Park, CA

15 91303-2944.

10

A method for generation of a 3D model of the dentition from an in-vivo scan of the

patient, and interactive computer-based treatment planning for orthodontic patients, is

described in published PCT patent application of OraMetrix, Inc., the assignee of this

invention, publication no. WO 01/80761, the contents of which are incorporated by reference

20 herein.

Other background references related to capturing three dimensional models of

dentition and associated craniofacial structures include S.M. Yamany and A.A. Farag, "A

System for Human Jaw Modeling Using Intra-Oral Images" in Proc. IEEE Eng. Med. Biol.

Soc. (EMBS) Conf., Vol. 20, Hong Kong, Oct. 1998, pp. 563-566; and M. Yamany, A.A.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

Farag, David Tasman, A.G. Farman, "A 3-D Reconstruction System for the Human Jaw

Using a Sequence of Optical Images," IEEE Transactions on Medical Imaging, Vol. 19, No.

5, May 2000, pp. 538-547. The contents of these references are incorporated by reference

herein.

5

10

15

20

The technical literature further includes a body of literature describing the creation of

3D models of faces from photographs, and computerized facial animation and morphable

modeling of faces. See, e.g., Pighin et al., Synthesizing Realistic Facial Expression from

Photographs, Computer Graphics Proceedings SIGGRAPH '98, pp. 78-94 (1998); Pighin et

al., Realistic Facial Animation Using Image-based 3D Morphing, Technical Report no.

UW-CSE-97-01-03, University of Washington (May 9, 1997); and Blantz et al., A Morphable

Model for The Synthesis of 3D Faces, Computer Graphics Proceedings SIGGRAPH '99

(August, 1999). The contents of these references are incorporated by reference herein.

The present invention is directed to an effective, computer-based, integrated and

interactive orthodontic treatment planning system that provides the necessary tools to allow

the orthodontist to quickly and efficiently design a treatment plan for a patient. The present

invention also provides a treatment planning system in which the orthodontist-derived

parameters for the treatment can be translated into a design of the treatment. The preferred

embodiment integrates 2D and 3D images to drive effective treatment planning. Intelligence

is built into the system whereby predefined therapeutic strategies, such as extraction,

interproximal reduction, distal movement of molars, can have associated value sets

predefined by the clinician that are used to drive the appropriate set-up automatically. Such

predefined therapeutic strategies could be entered via convenient user interface tools, such as

5

by templates.

The treatment design as described herein also allows for real-time communication of

the treatment plan to occur with the patient, or transmitted over a communications link and

shared with a colleague or remote appliance manufacturing facility. Alternatively, the

treatment planning can be performed remotely and a digital treatment plan sent to the

5 orthodontist for review, interactive modification, or approval.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

SUMMARY OF THE INVENTION

In a first aspect, an orthodontic treatment planning workstation is provided

comprising a computing platform having a graphical user interface, a processor and a

computer storage medium containing digitized records pertaining to a patient. The digitized

records include image data. The computer storage medium further includes a set of software

instructions providing graphical user interface tools for providing a user with access to the

digitized records for planning orthodontic treatment of a patient. The set of instructions

include:

5

10

15

20

a) treatment plan instructions providing graphical user interface tools for

allowing the user to interactively create a proposed set-up for treatment of the patient, the

proposed set-up comprising a proposed three-dimensional position of the teeth of the

upper and lower arches of the patient in a post-treatment condition; and

b) evaluation instructions providing a series of predetermined steps for guiding a

user to interactively evaluate the proposed set-up. The predetermined steps comprise steps

for 1) evaluation of the proposed set-up against boundary conditions for treatment of the

patient, the boundary conditions including a least a midline, an occlusal plane, a fixed

reference object (e.g. tooth or other anatomical structure that remains fixed) and an arch

form, and 2) evaluation of whether the tooth positions in both arches, and the inter-arch

relationship, of the proposed set-up are essentially ideal for treatment of the patient, that is,

correspond to the treatment goals of the patient. As used herein, the phrase "essentially ideal

for treatment of the patient" means that the proposed set up satisfies the practitioner's

objectives or goals for treatment of the patient, in other words, the proposed set up meets the

objectives for treatment of the patient, such as interarch relationship, tooth position, etc.,

midline placement, overbite and overjet, whatever they may be for a given patient. This

7

invention recognizes that there may be more than one possible "ideal" setup for a given

patient, as different practitioners may arrive at unique solutions for treatment for the patient,

and that for any given practitioner there may be more than one proposed setup that meets the

objectives for treatment of the patient, any one of which would satisfy the term "essentially

ideal for treatment of the patient." The term "ideal" is thus not used in an absolute sense of

the word.

10

15

20

In a related aspect, a computerized method of planning treatment for an orthodontic

patient is provided, comprising the step of providing an orthodontic treatment planning

workstation comprising a computing platform having a graphical user interface, a processor

and a computer storage medium containing digitized records pertaining to a patient. The

digitized records including image data. The workstation includes software instructions

providing graphical user interface tools for access to the digitized records and for planning

orthodontic treatment of a patient. The method continues with a step of using the workstation

to generate a proposed set-up for treating the patient, the proposed set-up comprising a

proposed three-dimensional position of the teeth of the upper and lower arches of the patient

in a post-treatment condition. The method further continues with a step of conducting an

evaluation of the proposed set-up, the evaluation prompted by computer instructions

providing a series of predetermined steps for guiding a user to interactively evaluate the

proposed set-up. The predetermined steps comprise steps for 1) evaluation of said proposed

set-up against boundary conditions for treatment of the patient, the boundary conditions

including a least a midline and an occlusal plane, and 2) evaluation of whether the tooth

positions in both arches, and the inter-arch relationship, of the proposed set-up are essentially

8

ideal for treatment of the patient.

BRIEF DESCRIPTION OF THE DRAWINGS

Presently preferred embodiments of the invention are described below in reference to

the appended drawings, wherein like reference numerals refer to like elements in the various

views, and in which:

5

10

15

Figure 1 is block diagram of a system for creating a three-dimensional virtual patient

model and for diagnosis and planning treatment of the patient.

Figure 2 is a flow chart showing a method of three-dimensional face creation from

scanning systems, which may be executed in software in the computer system of Figure 1.

Figure 3 is a flow chart showing an alternative method of three-dimensional face

model face creation using a plurality of possible input image or data formats, which may be

executed in software in the computer system of Figure 1.

Figure 4 is a flow chart showing a method of creating a complete textured three-

dimensional model of teeth.

Figures 4A-4E show a technique for combining 2D color photographs with 3D tooth

data to created textured (colored) 3D tooth models.

Figure 5 is a screen shot of the user interface of Figure 1 showing a three-dimensional

face model and a three-dimensional tooth model, in separate coordinate systems (i.e., prior to

registration or superposition of the two relative to each other). Figure 5 also shows a

plurality of icons, which, when activated, provide tools for manipulating the models shown in

20 the Figure.

Figure 6 is a screen shot showing one possible method of placement of the lower jaw

3D data into the face data coordinate system using corresponding points that are common to

9

each data set.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

Figure 7 is a screen shot showing the face data and the lower jaw 3D data in a

common coordinate system (the face coordinate system of Figures 5 and 6).

Figure 8 is a screen shot showing the face data and skull data obtained from a CT scan

in a common coordinate system.

Figure 9 is a screen shot showing face data and skull data superimposed on X-ray data

obtained from the patient.

Figure 10 is a screen shot showing the superposition of skull and face data with X-

Ray data.

5

10

20

Figures 11A-11E are a series of views of a digital model of an orthodontic patient

obtained, for example from CT scan, photographs, or intra-oral scanning with a hand-held 3D

scanner.

Figure 12 is a diagram illustrating a technique for scaling orthodontic data obtained

from an imaging device, such as a camera, to the actual anatomy of the patient.

Figure 13 is a diagram showing an alternative scaling method similar to that shown in

15 Figure 12.

Figure 14 is an illustration of an X-ray of a set of teeth and adjacent bone.

Figure 15 is an illustration of scaling the X-ray data of the tooth to the actual size of

the tooth to produce a scaled digital model of the tooth.

Figures 16A-16C are an illustration of a method of determining orientation reference

points in a digital model of a patient.

Figure 17 is an illustration of a method of mapping the orientation reference points of

Figures 16A-16C to a three-dimensional coordinate system.

Figure 18 is an illustration of a method of mapping the orientation reference points of

Figures 16A-16C to a three-dimensional coordinate system.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

Figure 19 is a more detailed block diagram of the treatment planning software

executed by the workstation of Figure 1.

Figure 20 is an illustration of the integration of the patient data acquisition, treatment

planning and appliance design functions that are facilitated by a preferred embodiment of the

5 unified workstation.

Figures 21-58 illustrate screen shots from the workstation of Figure 1, showing

various treatment planning features that are provided by the treatment planning software of

Figure 19 in a presently preferred embodiment.

Figure 59 is a diagram showing the types of evaluation that can be performed in this

invention.

10

15

Figures 60-91 illustrate various screen shots on the workstation of Figure 1 (or on a

remote workstation), showing various proposed set-up evaluation features that are provided

by the evaluation instructions executed in the workstation in a preferred embodiment.

Figure 92-94 are an index of icons that may be provided to the user on the display of

the workstation in various treatment planning and setup evaluation phases or modules in the

software, with the icons associated with the functions described in the index. In instances

where the different icons appear to be the same in these Figures, different colors are used for

lines or other features in the icons, and this variation in color is not reflected in the black and

white drawings used for Figures 92-94.

Figure 95 shows a screen shot of the user interface of the workstation, after the user

has selected individual brackets for use in treating the patient. The user has navigated to a

screen display in which the user is shown the particular brackets for each tooth. When they

highlight a bracket, a virtual model of the bracket is displayed and the prescription of the

bracket can be evaluated against a normative database.

25

20

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

DETAILED DESCRIPTION OF PREFERRED EMBODIMENTS

Before describing the treatment planning features of this invention in detail, an

overview of a unified workstation will be set forth initially. The workstation, in a preferred

embodiment, provides software features that create two dimensional and/or three-dimensional

virtual patient model on a computer, which can be used for purposes of treatment planning in

accordance with a presently preferred embodiment.

Many of the details and computer user interface tools which a practitioner may use in

adjusting tooth position, designing appliance shape and location, managing space between

teeth, and arriving at a finish tooth position using interaction with a computer storing and

displaying a virtual model of teeth are set forth in the prior application serial no. 09/834,412

filed April 13, 2001, and in published OraMetrix patent application WO 01/80761, the

contents of which are incorporated by reference herein. Other suites of tools and functions

are possible and within the scope of the invention. Such details will therefore be omitted

from the present discussion.

10

15

20

25

General Description

A unified workstation environment and computer system for diagnosis, treatment

planning and delivery of therapeutics, especially adapted for treatment of craniofacial

structures, is described below. In one possible example, the system is particularly useful in

diagnosis and planning treatment of an orthodontic patient. Persons skilled in the art will

understand that the invention, in its broader aspects, is applicable to other craniofacial

disorders or conditions requiring surgery, prosthodontic treatment, restorative treatment, etc.

A presently preferred embodiment is depicted in Figure 1. The overall system 100

includes a general-purpose computer system 10 having a processor (CPU 12) and a user

12

interface 14, including screen display 16, mouse 18 and keyboard 20. The system is useful

for planning treatment for a patient 34.

5

10

15

20

The system 100 includes a computer storage medium or memory 22 accessible to the

general-purpose computer system 10. The memory 22, such as a hard disk memory or

attached peripheral devices, stores two or more sets of digital data representing patient

craniofacial image information. These sets include at least a first set of digital data 24

representing patient craniofacial image information obtained from a first imaging device and

a second set of digital data 26 representing patient craniofacial image information obtained

from a second image device different from the first image device. The first and second sets

of data represent, at least in part, common craniofacial anatomical structures of the patient.

At least one of the first and second sets of digital data normally would include data

representing the external visual appearance or surface configuration of the face of the patient.

In a representative and non-limiting example of the data sets, the first data set 24

could be a set of two dimensional color photographs of the face and head of the patient

obtained via a color digital camera 28, and the second data set is three-dimensional image

information of the patient's teeth, acquired via a suitable scanner 30, such as a hand-held

optical 3D scanner, or other type of scanner. The memory 22 may also store other sets 27 of

digital image data, including digitized X-rays, MRI or ultrasound images, CT scanner etc.,

from other imaging devices 36. The other imaging devices need not be located at the location

or site of the workstation system 100. Rather, the imaging of the patient 34 with one or other

imaging devices 36 could be performed in a remotely located clinic or hospital, in which case

the image data is obtained by the workstation 100 over the Internet 37 or some other

communications medium, and stored in the memory 22.

The system 100 further includes a set of computer instructions stored on a machine-

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

readable storage medium. The instructions may be stored in the memory 22 accessible to the

general-purpose computer system 10. The machine-readable medium storing the instructions

may alternatively be a hard disk memory 32 for the computer system 10, external memory

devices, or may be resident on a file server on a network connected to the computer system.

the details of which are not important. The set of instructions, described in more detail

below, comprise instructions for causing the general computer system 10 to perform several

functions related to the generation and use of the virtual patient model in diagnostics,

therapeutics and treatment planning.

10

15

20

These functions include a function of automatically, and/or with the aid of operator

interaction via the user interface 14, superimposing the first set 24 of digital data and the

second set 26 of digital data so as to provide a composite, combined digital representation of

the craniofacial anatomical structures in a common coordinate system. This composite,

combined digital representation is referred to herein occasionally as the "virtual patient

model," shown on the display 16 of Figure 1 as a digital model of the patient 34. Preferably,

one of the sets 24, 26 of data includes photographic image data of the patient's face, teeth and

head, obtained with the color digital camera 28. The other set of data could be intra-oral 3D

scan data obtained from the hand-held scanner 30, CT scan data, X-Ray data, MRI, etc. In

the example of Figure 1, the hand-held scanner 30 acquires a series of images containing 3D

information and this information is used to generate a 3D model in the scanning node 31, in

accordance with the teachings of the published PCT application of OraMetrix, PCT

publication no. WO 01/80761, the content of which is incorporated by reference herein.

Additional data sets are possible, and may be preferred in most embodiments. For example

the virtual patient model could be created by a superposition of the following data sets: intra-

oral scan of the patient's teeth, gums, and associated tissues, X-Ray, CT scan, intra-oral color

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

photographs of the teeth to add true color (texture) to the 3D teeth models, and color

photographs of the face, that are combined in the computer to form a 3D morphable face

model. These data sets are superimposed with each other, with appropriate scaling as

necessary to place them in registry with each other and at the same scale. The resulting

representation can be stored as a 3D point cloud representing not only the surface on the

patient but also interior structures, such as tooth roots, bone, and other structures. In one

possible embodiment, the hand-held in-vivo scanning device is used which also incorporates

a color CCD camera to capture either individual images or video.

The software instructions further includes a set of functions or routines that cause the

user interface 16 to display the composite, combined digital three-dimensional representation

of craniofacial anatomical structures to a user of the system. In a representative embodiment,

computer-aided design (CAD)-type software tools are used to display the model to the user

and provide the user with tools for viewing and studying the model. Preferably, the model is

capable of being viewed in any orientation. Tools are provided for showing slices or sections

through the model at arbitrary, user defined planes. Alternatively, the composite digital

representation may be printed out on a printer or otherwise provided to the user in a visual

form.

10

15

20

The software instructions further include instructions that, when executed, provide the

user with tools on the user interface 14 for visually studying, on the user interface, the

interaction of the craniofacial anatomical structures and their relationship to the external,

visual appearance of the patient. For example, the tools include tools for simulating changes

in the anatomical position or shape of the craniofacial anatomical structures, e.g., teeth, jaw,

bone or soft tissue structure, and their effect on the external, visual appearance of the patient.

The preferred aspects of the software tools include tools for manipulating various parameters

15

such as the age of the patient; the position, orientation, color and texture of the teeth;

reflectivity and ambient conditions of light and its effect on visual appearance.

elements of the craniofacial and dental complex can be analyzed quickly in either static

format (i.e., no movement of the anatomical structures relative to each other) or in an

dynamic format (i.e., during movement of anatomical structures relative to each other, such

as chewing, occlusion, growth, etc.). To facilitate such modeling and simulations, teeth

may be modeled as independent, individually moveable three dimensional virtual objects,

using the techniques described in the above-referenced OraMetrix published PCT application,

WO 01/80761.

10

15

The workstation environment provided by this invention provides a powerful system

and for purposes of diagnosis, treatment planning and delivery of therapeutics. For example,

the effect of jaw and skull movement on the patient's face and smile can be studied.

Similarly, the model can be manipulated to arrive at the patient's desired feature and smile.

From this model, and more particularly, from the location and position of individual

anatomical structures (e.g., individual tooth positions and orientation, shape of arch and

position of upper and lower arches relative to each other), it is possible to automatically back

solve for or derive the jaw, tooth, bone and/or soft tissue corrections that must be applied to

the patient's initial position, which might be pre-treatment position or position at any other

time during treatment, to provide the desired result. This leads directly to a patient treatment

20 plan.

These simulation tools, in a preferred embodiment, comprise user-friendly and

intuitive icons 35 that are activated by a mouse or keyboard on the user interface of the

computer system 10. When these icons are activated, the software instruction provide pop-

up, menu, or other types screens that enable a user to navigate through particular tasks to

16

highlight and select individual anatomical features, change their positions relative to other

structures, and simulate movement of the jaws (chewing or occlusion). Examples of the

types of navigational tools, icons and treatment planning tools for a computer user interface

that may be useful in this process and provide a point of departure for further types of

displays useful in this invention are described in the patent application of Rudger Rubbert et

al., Serial No. 09/835,039 filed April 13, 2001, the contents of which are incorporated by

reference herein.

10

15

20

The virtual patient model, or some portion thereof, such as data describing a three-

dimensional model of the teeth in initial and target or treatment positions, is useful

information for generating customized orthodontic appliances for treatment of the patient.

The position of the teeth in the initial and desired positions can be used to generate a set of

customized brackets, and customized flat planar archwire, and customized bracket placement

jigs, as described in the above-referenced Andreiko et al. patents. Alternatively, the initial

and final tooth positions can be used to derive data sets representing intermediate tooth

positions, which are used to fabricate transparent aligning shells for moving teeth to the final

position, as described in the above-referenced Chisti et al. patents. The data can also be used

to place brackets and design a customized archwire as described in the previously cited

application serial no. 09/835,039.

To facilitate sharing of the virtual patient model among specialists and device

manufacturers, the system 100 includes software routines and appropriate hardware devices

for transmitting the virtual patient model or some subset thereof, and a proposed set-up for

treatment of the patient, over a computer network. The treatment plan developed in the

workstation could be either evaluated locally, on the workstation itself, or else at a remote

workstation (such as at an appliance manufacture site of practitioners that perform

17

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive

Chicago, IL 60606

(312) 913-0001

evaluations for a service). In this latter situation, the same interactive treatment planning

software is installed at the remote site, along with the evaluation instructions for evaluation of

the set-up, as described herein.

The workstation's software instructions are preferably integrated with a patient

management program having a scheduling feature for scheduling appointments for the

patient. The patient management program provides a flexible scheduling of patient

appointments based on progress of treatment of the craniofacial anatomical structures. The

progress of treatment can be quantified. The progress of treatment can be monitored by

periodically obtaining updated three-dimensional information regarding the progress of

treatment of the craniofacial features of the patient, such as by obtaining updated scans of the

patient and comparison of the resulting 3D model with the original 3D model of the patient

prior to initiation of treatment.

10

15

20

Thus, it is contemplated that system described herein provides a set of tools and data

acquisition and processing subsystems that together provides a flexible, open platform or

portal to a variety of possible therapies and treatment modalities, depending on the preference

of the patient and the practitioner. For example, a practitioner viewing the model and using

the treatment planning tools may determine that a patient may benefit from a combination of

customized orthodontic brackets and wires and removable aligning devices. Data from the

virtual patient models is provided to diverse manufacturers for coordinated preparation of

customized appliances. Moreover, the virtual patient model and powerful tools described

herein provide a means by which the complete picture of the patient can be shared with other

specialists (e.g., dentists, maxilla-facial or oral surgeons, cosmetic surgeons, other

orthodontists) greatly enhancing the ability of diverse specialists to coordinate and apply a

diverse range of treatments to achieve a desired outcome for the patient. In particular, the

18

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

overlay or superposition of a variety of image information, including 2D X-Ray, 3D teeth

image data, photographic data, CT scan data, and other data, and the ability to toggle back

and forth between these views and simulate changes in position or shape of craniofacial

structures, and the ability to share this virtual patient model across existing computer

networks to other specialists and device manufacturers, allows the entire treatment of the

patient to be simulated and modeled in a computer. Furthermore, the expected results can be

displayed before hand to the patient and changes made depending on the patient input.

With the above general description in mind, additional details of presently preferred

components and aspects of the inventive system and the software modules providing the

functions referenced above will be described next.

Capture of image information

10

15

20

25

The creation of the virtual patient model uses the capture and storage of at least two

different digital sets of image data of the patient. The image sets will typically represent, at

least in part, overlapping craniofacial anatomical structures so that a superposition of them in

a common three-dimensional coordinate system may occur.

The type of image data that will be obtained will vary depending on the available

image acquisition devices available to the practitioner. Preferably, the system employs

software simulation of changes in shape or position of craniofacial structures (e.g., teeth or

jaw) on the visual appearance, e.g., smile, of the patient. Accordingly, at least one of the data

sets will normally include data regarding the surface configuration of the face and head. A

commercially available digital CCD camera 28 (Fig. 1), e.g., camera available from Sony or

Canon, can be used to obtain this information. Preferably, the image data is color image data.

The data sets are obtained by photographing the patient's head and face at various viewing

19

angles with the camera and storing the resulting image files in the memory of the computer.

These images can provide a basis for creating a morphable face model.

The image data regarding the patient's exterior appearance can be obtained through

other means including via scanning of the head and face of the patient via the hand-held 3D-

scanner 30 described in the published OraMetrix PCT application, referenced previously. If

this approach is used, it may be beneficial to apply a thin layer of non-toxic, opaque and

reflective substance to the skin prior to scanning to insure adequate data capture by the hand-

held scanner. A suitable opaquing substance is described in the patent application of Nancy

Butcher et al. serial no. 10/099,042 filed March 14, 2002, entitled "Method for Wet-Field

Scanning," the contents of which are incorporated by reference herein. In operation, the

scanner captures a sequence of overlapping images of the surface of the patient as the scanner

is held by the hand and moved about the face. The set of images can be obtained in only a

few minutes. Each image is converted to a set of X, Y and Z coordinate positions comprising

a cloud of points representing the surface of the face. The point clouds from each image are

registered to each other to find a best fit to the data. The resulting registered point cloud is

then stored in the memory as a virtual three-dimensional object. The construction, calibration

and operation of the scanner, and the manner of converting scanned data to point clouds and

registering three-dimensional point clouds to form a three-dimensional object is described at

length in the published PCT application of OraMetrix and therefore omitted from the present

discussion for the sake of brevity. Other types of scanners or coordinate measuring

instruments could be used in less preferred embodiments, such as the scanning devices in the

Yamany et al. articles referenced previously.

Aside from surface data of the patient obtained by the camera 28 or 3D scanner 30,

the system typically will include the capture of additional data representing the teeth of the

20

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

5

10

15

patient, and also capture of additional data representing craniofacial structures not visible to

the naked eye using other imaging devices 36 (Figure 1). For example, the system will

acquire digitized images from an X-ray machine capturing X-ray photographs of the patient's

head, jaw, teeth, roots of teeth, and other craniofacial structures. These photographs are

digitized and stored in the memory of the computer system.

As other possible examples, three-dimensional magnetic resonance images of the

patient's head or jaws are obtained and stored in the memory. Other examples include

images acquired from a computed tomography (CT) scanner, ultrasound imager, or other type

of imaging device.

5

10

15

20

While the above discussion has described how 3D image of the face can be obtained

from a three-dimensional scanner, there are other possibilities that may be used in the

practice of alternative embodiments. One such alternative is creating a 3D virtual face from a

series of 2-D color photographs. This technique is known and described in Pighin et al.,

Synthesizing Realistic Facial Expression from Photographs, Computer Graphics Proceedings

SIGGRAPH '98, pp. 78-94 (1998); Pighin et al., Realistic Facial Animation Using Image-

based 3D Morphing, Technical Report no. UW-CSE-97-01-03, University of Washington

(May 9, 1997); and Blantz et al., A Morphable Model for The Synthesis of 3D Faces,

Computer Graphics Proceedings SIGGRAPH '99 (August, 1999), the contents of which are

incorporated by reference herein. Basically, in this alternative, two-dimensional color

pictures of the face are taken which are converted automatically to a textured 3 dimensional

model using a 'morphable model' technique. Here, the phrase "textured 3 dimensional

model" is used in the particular sense of a colorized three-dimensional object, with the word

21

"texture" synonymous with color data, as that term is used in this particular art.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

Morphable models can be built based on various known approaches such as optic

flow algorithms or active model matching strategy, or a combination of both. One approach

is to scan a set of 2D faces. A shape vector containing 3D vertices and texture vector

containing RGB color values of each vertex represents the geometry of the face. Each face is

divided into sub regions such as eyes, nose, mouth etc. Blending the sub-regions at the

borders generates the complete 3D face. Automatic matching and estimating 3D face of a 2D

color image from morphable model is carried out as follows:

New Shape (Sn) and texture (Tn) are computed as follows:

(1) $\operatorname{Sn} = \operatorname{Sa} + \sum \alpha s$;

5

10

15

20

(2) Tn = Ta + $\Sigma \beta t$, where Sa and Ta are the averages of Shape S and Texture T over all

the 3D face datasets; s & t are the eigenvectors of the covariance matrices; α and β are the

coefficients of the facial shape and texture for all the faces, and n is a sub-region index.

Rendering parameters p contain camera position, object scale, image plane rotation

and translation and light intensity. From Bayes decision theory, the set of parameters,

 (α,β,ρ) are determined with maximum posterior probability for getting a corresponding 3D

face from a 2D image.

Three-dimensional image data sets of the upper and lower arches including upper and

lower teeth are preferably created with a 3D optical scanner 30, such as the OraMetrix hand-

held in-vivo scanner. If the 3D jaw model has no texture model, i.e., no color data, the

texture data can be extracted from the 2 dimensional colored picture of the upper and lower

jaw and mapped to the 3D coordinates on the jaw model using a cylindrical projection

technique. In this technique, a map is constructed in texture space, that for each point (u, v),

specifies a triangle whose cylindrical projection covers that point. The 3D point p

22

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

corresponding to point (u, v) in texture space is computed by intersecting a ray with the surface of the corresponding point in the 2D colored image.

Superposition or registration of the data sets

5

10

15

20

After the images of the face, craniofacial structures, X-rays, teeth etc. are obtained and stored in memory in digital form they are superimposed on each other (i.e., registered to each other via software in the workstation) to create a complete virtual patient model on the workstation. The superposition of the sets of image data may be developed as an automatic software process, or one in which there is user involvement to aid in the process. possible example, the three-dimensional textured model of the face is properly aligned with the 3D jaw model obtained from the intra-oral scan, 3D skull data from CT scan, and 2 dimensional X-rays to create a virtual patient model. For correct alignment of the data sets to each other, a preferred method executed by the software selects three or more corresponding points on the 3D jaw and the 3D face, and then computes a transformation matrix to re-orient the 3D face relative to the 3D jaw. This transformation matrix will contain the information needed to rotate and translate the 3D face relative to the 3D jaw in a best-fit manner to align the two to each other. Methods of calculation of transformation matrices to achieve registration are taught in the published PCT patent application of OraMetrix, Inc., cited previously. Similar methods are used for registering the CT scan data and X-ray data to the combined 3D face and jaw model. Once the superposition is achieved, the resulting model is displayed on the workstation user interface. The user is provided with tools for simulating movement or repositioning of craniofacial structures of the virtual patient, and the computer animates such movement or repositioning and shows the effect of such movement or repositioning on the external visual appearance of the patient.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

An example of registering scan data of a jaw from an intra-oral scan to a 3D face model using human interaction is shown in Figures 2-7. Figure 2 is a flow chart showing a software method of three-dimensional face creation from scanning systems, which may be executed in software in the computer system 10 of Figure 1. There are two possible approaches for creating the 3D face, one using a color digital camera 28 (Figure 1) and another using scanning of the face using the hand held scanner 30 and scanning node 31 In the situation in which a color digital camera is used, at step 40 a set 24 of 2D digital color photographic images of the face and head are obtained and stored in the memory 22 of Figure 1. The set 24 of images is supplied to a module 42 which creates a virtual 3D face using an active model matching strategy, using the techniques known in the art and described in Pighin et al., Synthesizing Realistic Facial Expression from Photographs, Computer Graphics Proceedings SIGGRAPH '98, pp. 78-94 (1998); Pighin et al., Realistic Facial Animation Using Image-based 3D Morphing, Technical Report no. UW-CSE-97-01-03, University of Washington (May 9, 1997); and Blantz et al., A Morphable Model for The Synthesis of 3D Faces, Computer Graphics Proceedings SIGGRAPH '99 (August, 1999). This 3D face model is then stored in the hard disk memory of the computer 10, as indicated at process module 44.

In alternative embodiments, a 3D scanning of the face using a laser or 3D optical scanner is performed, as indicated at 46. The 3D model is provided to a module 48 which creates a morphable model of the face and head with an optic flow algorithm. This morphable model is provided to the module 42 for creating a 3D face. At step 50, the software inquires as to whether a morphable 3D face is available, and if not the processing of module 42 executes, in which a 3D morphable model of the face is created. If a morphable 3D face is already available, the software inquires at step 54 as to whether texture (color)

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

10

15

information is available to add to the 3D face. (Note that in many 3D scanner systems there

is no acquisition of color information, only spatial information). If color information is not

available, the processing proceeds to module 56. In module 56, the color data is provided to

the 3D model to create a 3D color morphable virtual model. The color data is supplied from

the digital photographs of the patient, obtained at step 40. The texture information is

supplied to the 3D model from the scanner using a cylindrical projection technique in module

56 (or by using any other known technique). The textured, morphable 3D model of the face

and head is stored as indicated at module 44.

An alternative software method or process for creating a 3D model of the face is

shown in Figure 3. The method involves the acquisition of a 3D color face model 52 (using

for example the techniques of Figure 2), the acquisition of 3D color model of the teeth 54,

and the acquisition of a model 56 of the skull using a CT scanner. These three models are

supplied to a module 58 which performs an aligning transformation on the data sets from

each of these modules. The aligning transformation process 58 basically scales and provides

the necessary X, Y and Z translations and rotations to place the data sets into a common

coordinate system such that common anatomical structures overlap each other. The

complete 3D face model is stored as indicated at step 60 and then supplied to an overlay

transformation module 66. The overlay transformation module 66 obtains a set of 2D color

face photographs 62 and X-Rays 64, and overlays them to the complete 3D face model to

result in a combined, composite model of the face, skull, teeth, and associated tooth roots,

bone and other anatomical data. This composite representation of the patient is stored in a

database 68 for the system 100.

Figure 4 shows a process that can be used to combine 3D scan data with 2D color

photographs to create a 3D color model of the teeth. In step 70, the teeth are scanned with

25

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

10

15

the hand-held intra-oral scanner 30 of Figure 1. The resulting data represent a 3D model of

the dentition, which is stored in the computer 10. This process is described in the published

PCT application of OraMetrix, Inc. cited previously. At step 72, 2D color photographs of

the teeth are obtained with a color digital camera. In one possible embodiment, the hand-

held scanner 30 includes a video camera that obtains a continuous stream of color video

frames separate and apart from the acquisition of 3D image data. The color photographs of

the dentition at step 72 could be obtained in this manner.

At step 74, a 3D textured model of the teeth is created using a cylindrical projection

technique. Basically, in this technique, the color data from the color photographs is projected

onto the tooth data. The tooth data can be represented as triangular surfaces, with the vertices

of each triangle being adjacent points in a point cloud defining the surface of the tooth. The

color is projected on the surfaces, and each surface is assigned a value associated with a

particular color. The result is a 3D color model of the teeth at step 76.

Figures 4A-4E show several screen displays from a user interface of the unified

workstation that illustrate the process of texture mapping a 3D object (here, teeth) by

projection of color data from a 2D photograph. After a patient's dentition is scanned, the

virtual teeth and gingiva for both upper and lower arches are represented as a single surface,

in the present example a triangle mesh surface. Figure 4A shows a 2D digital photograph of

teeth/gingivae 71 displayed in a graphical window 73 along with a 3D virtual model of the

teeth 75 to one side. The 2D digital photograph 71 is scaled up or down in size as

necessary so as to be approximately the same in scale (size) as the 3D model of the teeth 75.

This is accomplished using any suitable icons or mouse action, such as clicking on the 2D

photograph and scrolling up or down with the mouse to change the size of the 2D image so

that it matches the size of the 3D model. Figure 4B shows the surface of the teeth and

26

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

10

15

gingivae of the 3D virtual model 75 in greater detail. The surface of the model 75 comprises a

set of minute interconnecting triangle surfaces, with the vertices of the triangle surfaces being

points that represent the surface. This is only one possible format for representing the

surface of a 3D object.

5

10

15

20

After the 2D photograph and 3D model have been scaled, a translation is performed

so as to overlap the 3D model and the 2D photograph. Figure 4C shows the 2D picture 71

transformed by scaling and translation such that it is superimposed on the 3D model 75. This

superposition could be performed manually or automatically. For example, the user can

click and drag the 2D digital photograph 71 and manually move it using the mouse so that it

overlaps exactly the 3D model 75. The color information in the 2D photograph 71 is

projected and mapped to the individual triangle surfaces forming the lower jaw and upper jaw

of the 3D model 75 using, for example, a projection algorithm. The result, a textured 3D

model, is shown in Figure 4D. Figure 4E shows textured 3D model after rotation on the user

interface.

Occlusal and lingual 2-D color photographs of each jaw are also obtained and texture

data is mapped to the surface data. The result is a complete true color 3D model of the teeth

of both arches.

Figure 5 is an illustration of a screen display on the user interface of the computer 10.

The display shows a 3D morphable model 102 of patient on the left hand side of the display,

in a given arbitrary coordinate system X, Y, Z. The morphable model 102 is obtained, for

example, from color photographs using the techniques described previously. A three-

dimensional model 104 of teeth of the patient is shown on the right hand side of the screen.

The 3D model of the teeth 104 can be obtained from intra-oral scanning using the scanner 30

of Figure 1, from a laser scan of a physical model of the dentition obtained from an

27

impression, from a coordinate measuring device or some other source. The source is not

particularly important. The 3D model of the teeth 104 is shown in a separate coordinate

system X', Y', Z'. Screen display includes various icons 35 that allow the user to position

the tooth model 104 relative to the morphable model 102 in order to combine the two in a

common coordinate system and construct a composite model.

In Figure 6, the user has activated an "Align References" icon, which causes the

screen display to show the box 106 on the left hand side of the screen. The user is provided

with the option to pick points that represent anatomical structures that are common to both

the morphable model 102 and the 3D tooth model 104. In this particular situation, the user

has selected with the mouse two points on the lower arches which lie at the intersection of the

teeth and the gums. These two points are shown as a triangle 108 and a square 110.

Obviously, other points could be chosen. The user then clicks on the "Apply" tab 112. The

result is shown in Figure 7, in which the 3D tooth model 104 is combined with the morphable

face 102 model to produce a combined virtual patient model 34.

In the example of Figures 5-7, the morphable model 102 was already scaled to the

same scale as the tooth model 104. In other words, the data representing the morphable face

model indicates that the spatial dimensions of the teeth in the morphable face model is

substantially the same as the spatial dimensions of the virtual tooth model 104. Methods of

performing scaling are described below.

Figure 8 is an illustration of an alternative embodiment of a virtual patient model 34.

In this embodiment, the model 34 is a combination of data 102 representing a morphable

face, obtained from a plurality of 2D color photographs, and skull data 114 obtained from a

The two sets of data are shown in a common coordinate system, appropriately

scaled. The manner of combining the two data sets can be using the approach described in

28

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

10

15

Figures 6 and 7. Alternatively, the user could click and drag using a mouse one or the other of

the data sets and manually position it until it is in the correct position. As yet another

alternative, the software could be programmed to find common overlapping features (such as

for example teeth) using surface matching algorithms, and then position the CT scan model

relative to the face model such that the common features overlap exactly.

Figure 9 is a screen shot of yet another possible embodiment of a virtual patient

model. This model combines face data 102 from a morphable face model (obtained from 2D

color photographs), skull data 114 from a CT scan, and X-ray data 116 obtained from a set of

digital X-Rays of the patient. The manner of creating the virtual patient model can be for

example using the procedure of Figure 3 and Figures 6-7. The morphable face model is

aligned relative to the CT scan data either automatically or using some human interaction.

The 2D X-Ray data can be morphed into 3D digital data using the morphable model

algorithms cited previously. Alternatively, the 2D X-Ray data can be combined with 3D

optical scan data of crowns of the teeth to create a combined X-Ray/3D tooth model, which is

then combined with the CT/morphable face model. This process may be optimized by using

virtual template tooth roots, which are modified to fit the X-ray data, and then this combined

3D root model is combined with crown data to produce a complete set of 3D teeth, including

roots. This combined model is merged into the CT scan/morphable face model using the

techniques of Figures 6 and 7 (selecting common points then using the "Apply Points" icon,

Figure 7, item 112), using click and drag techniques, or any other appropriate registration or

transformation technique.

Once the virtual model is created, the user is provided with tools that allow the user to

hide one or more image data in order to study certain features. Furthermore, the user is

provided with navigation tools with the ability to manipulate the model so as to view it from

29

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

10

15

any user-specified perspective. For example, in Figure 10 a screen shot is shown of the

superposition of skull data 114 with X-Ray data 116. In this example, complete 3D models

of the teeth 116 are created from X-ray data using the algorithms described previously.

Alternatively, the complete 3D tooth models 116 are created from combining X-Ray data

with 3D models of teeth obtained by a scan of the crowns of the teeth (using the scanner 30

of Figure 1 or from a laser scan of a physical model of the dentition), and/or with the use of

template tooth roots that are modified to match the X-ray data.

Scaling of Data

10

15

20

When digital image data from multiple sources are combined or superimposed

relative to each other to create a composite model, it may be necessary to scale data from one

set to the other in order to create a single composite model in a single coordinate system in

which the anatomical data from both sets have the same dimensions in three-dimensional

space. Hence, some scaling may be required. This section describes some approaches to

scaling that may be performed in one possible embodiment of the invention.

Figure 11A-11E are views of scan data 200 representing the front, right side, left side,

lower and upper arches of a patient. The data includes the teeth 202 and the gingiva 204.

Figure 12 illustrates a technique of scaling the orthodontic data to match the actual

orthodontic size. Depending on of the scanning technique, the orthodontic data will not

completely reproduce the exact size of the teeth and other portions of the orthodontic

structure. To facilitate the accuracy of the three-dimensional digital model, at least one tooth

206 can be marked utilizing one or more markings 208. The marking is done prior to

obtaining the orthodontic data. Once the orthodontic data for the tooth 206 is obtained,

scaling reference points 210 are also obtained. The scaling reference points are the points in

the scan data that represent the image of the markings 208. A determination of the

30

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

differences between the scaling reference points 210 and the actual markings 208 determine a

scaling factor 212. As one of average skill in the art will readily appreciate, having the actual

markings 208 and the scaling reference points 210, a variety of mathematical operations may

be used to determine the scaling factor 212. For example, the coordinate differences

(distance) between each of the vertices of the triangle may be utilized. As one of average

skill in the art will further appreciate, a different number of markings 208 may be utilized.

For example, two markings may be used or four markings may be used, etc. In addition,

more than one tooth may be marked with similar markings 208. Note that the markings may

be on the exterior of the patient, and a local triangulation technique may be used to obtain the

scaling factor. Further note that the scaling factor 212 determination is based on an

assumption that the scan data will have a linear error term in each of the x, y and z axis, such

that a single scaling factor is determined and used to scale each of the teeth as well as the

other aspects of the orthodontic structure of the patient.

Figure 13 illustrates an alternate marking technique for determining a scaling factor

for the orthodontic data. As shown, an actual tooth 206 is marked with a marking 208. The

marking 208 is of a substantial size so as to be adequately measured. Once the orthodontic

data is obtained, the orthodontic data of the tooth 202 and a corresponding scaling reference

point (area) 210 are used to determine the scaling factor 212. As one of average skill in the

art will readily appreciate, a simple mathematical function may be used to determine the

scaling factor 212 based on the size (diameter) difference between the actual marking 208

and the scaling reference point 210. As an alternative to marking as described with reference

to Figures 12 and 13, the actual tooth size, and the size of the model of the tooth, may be

measured and used to determine the scaling factor. Accordingly, the difference between the

31

actual tooth size and the size of the tooth in the scan data will constitute the scaling factor.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

10

15

When three-dimensional scanning of the type described in the published PCT

application of OraMetrix is used, scaling of the three-dimensional data is not needed as a

true, accurate and to scale three-dimensional image is obtained through the use of

triangulation. Likewise, a true three-dimensional image can be obtained using techniques

such as computed tomography. However, for video or photographic data, and for X-ray data,

scaling such as shown in Figures 12 and 13 may be needed in order to merge that data to

other data such as 3D scan data.

Figure 14 illustrates a two-dimensional representation of image data, such as a

graphical diagram of a radiographic image, such as an x-ray of a few teeth. In another

embodiment, the radiographic image can be a computed tomographic image volume. As

previously mentioned, the orthodontic data contains three-dimensional images of the surface

of the orthodontic structure. X-rays provide a more detailed view of the teeth and

surrounding hard and soft tissue as two dimensional image data. As shown in Figure 14, each

tooth 206 includes a crown 220, and a root 222 and is embedded in bone 224. Accordingly,

the orthodontic data 200 of Figure 11 only illustrates the crown 206 of the teeth. As such, a

complete three-dimensional model of the orthodontic patient requires the roots and bone to be

included.

10

15

20

Figure 15 illustrates a technique of using the scaled digital model 226 of the tooth's

crown to produce an integrated or composite digital model 228 of the tooth. In this

embodiment, the x-rayed data 230 of the tooth is used in comparison with the scaled digital

model 226 to determine a per tooth scaling factor. The scaled digital model 226 of the tooth

is positioned to be planar with the x-ray of the tooth 230. Having obtained the proper

orientation between the two objects, the per tooth scaling factor is determined and

subsequently used to generate the composite scaled digital model 228 of the tooth. In a

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive

(312) 913-0001

specific embodiment, the per tooth scaling factor is required for current x-ray technologies,

since x-rays produce a varying amount of distortion from tooth to tooth depending on the

distance of the tooth from the film, the angle of x-ray transmission, etc.

To more accurately map the two-dimensional images of a tooth onto the three-

dimensional model, multiple angles of the tooth should be used. Accordingly, a side, a front,

and a bottom view of the tooth should be taken and mapped to the scaled digital model of the

tooth. Note that the bone and other portions of the orthodontic structure are scaled in a

similar manner. Further note that MRI images, and any other images obtained of the

orthodontic patient, may also be scaled in a similar manner. A more complete

representation of the tooth roots may be obtained using standardized, template 3D virtual

tooth roots, applying the X-Ray data to the template tooth roots and modifying their shape

accordingly, and them applying the modified template tooth root to the scan data of the crown

to create a scaled, complete virtual tooth object including tooth roots.

Figures 16A-16C illustrate a graphical diagram of selecting orientation reference

points based on physical attributes of the orthodontic structure. The orientation reference

points 262 and 266 will be subsequently used to map the digital image of the orthodontic

structure into a three-dimensional coordinate system that will not change during the course of

treatment. In this example, the frenum 264 has been selected to be one of the orientation

reference points 266 and the rugae 260 has been selected as the other reference point 262.

The frenum 264 is a fixed point in the orthodontic patient that will not change, or change

minimally, during the course of treatment. As shown, the frenum 264 is a triangular shaped

tissue in the upper-portion of the gum of the upper-arch. The rugae 260 is a cavity in the roof

of the mouth 268 in the upper-arch. The rugae will also not change its physical position

through treatment. As such, the frenum 264 and the rugae 260 are fixed physical points in

33

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

10

15

the orthodontic patient that will not change during treatment. As such, by utilizing these as

the orientation reference points 262 and 266, a three-dimensional coordinate system may be

mapped thereto. Note that other physical attributes of the orthodontic patient may be used as

the orientation reference points 262 and 266. However, such physical points need to remain

constant throughout treatment. Accordingly, alternate physical points include the incisive

papilla, cupid's bow, the inter-pupillar midpoint, inter-comissural midpoint (e.g., between the

lips), inter-alar midpoint (e.g., between the sides of the nose), the prone nasale (e.g., the tip of

the nose), sub-nasale (e.g., junction of the nose and the lip), a dental mid-line point, a point

on the bone, a fixed bone marker such as an implant (e.g., a screw from a root canal, oral

surgery).

5

10

15

20

The x, y, z coordinate system may be mapped to the physical points on the digital

model of the orthodontic structure in a variety of ways. In one example, the origin of the x,

y, z coordinate system may be placed at the frenum 264, the z-axis aligned with reference to

the frenum and the rugae 260, and the x-axis is aligned with the midline of the upper and/or

lower arch. This is further illustrated in Figures 17 and 18. Note that an external positioning

system may be used to obtain the orientation reference points. For example, the patient may

sit in a chair at a specific location of an examination room that includes a triangulation

positioning system therein. As such, when the patient is scanned, the scanned images may be

referenced with respect to the room's triangulation positioning system.

Figure 17 illustrates a graphical representation of mapping the orientation reference

points 262 and 266 to the x-z plane of the three-dimensional x, y, z coordinate system. In this

illustration, orientation point 266, which corresponds to the frenum 264, is selected as the

origin of the x, y, z coordinate system. Note that any location may be selected as the origin

72. The orientation points 262 and 266 are used to determine an x, z plane orientation angle

34

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

262. Typically, the x, y, z coordinate system will be selected such that when looking at the

patient from a frontal view, the x direction will be to right of the patient, the y direction

towards the top of the patient's head and the z direction will be out away from the patient.

As one of average skill in the art will appreciate, the orientation of the x, y, z plane may be in

any orientation with respect to the reference points 262 and 266.

The x-y plane is mapped to the orientation reference point 262 and 266 as shown in

Figure 18. The orientation reference point 262 and 266 are used to generate an x-y plane

orientation angle 284. Based on the x-y plane orientation angle 284 and the x-z plane

orientation angle 262, a digital model of a tooth 270 may be positioned in three-dimensional

space with respect to the x, y, z coordinate system. As shown in Figures 17 and 18, the

digital model of the tooth 270 includes a tooth depth 278, an angle of rotation 276 with

respect to the x-z axis, an angle of rotation 282 with respect to the x-y plane, a positioning

vector 274 which is in a three-dimensional space, the length of the tooth including the crown

dimension, and the root dimension. Accordingly, each tooth is then mapped into the x, y, z

coordinate system based on the tooth's center, or any other point of the tooth, and the

dimensions of the digital model of the corresponding tooth. Once each tooth has been placed

into the x, y, z coordinate system, the digital model of the tooth is complete. Note that the

lower-arch is also referenced to the x, y, z coordinate system wherein the determination is

made based on the occlusal plane of the patient's orthodontic structure. Alternatively, the

lower-arch may include a separate three-dimensional coordinate system that is mapped to the

coordinate system of the upper-arch. In this latter example, fixed points within the lower-

arch would need to be determined to produce the lower arch's three-dimensional coordinate

35

system.

5

10

15

20

Treatment Planning

The computer or workstation 10 (Figure 1) that includes the software for generating

the patient model preferably includes interactive treatment planning software that allows the

user to simulate various possible treatments for the patient on the workstation and visualize

the results of proposed treatments on the user interface by seeing their effect on the visual

appearance of the patient, especially their smile. The interactive treatment planning

preferably provides suitable tools and icons that allow the user to vary parameters affecting

the patient. Such parameters would include parameters that can be changed so as to simulate

change in the age of the patient, and parameters that allow the user to adjust the color,

texture, position and orientation of the teeth, individually and as a group. The user

manipulates the tools for these parameters and thereby generates various virtual patient

models with different features and smiles. The patient models are displayed on the user

interface of the workstation where they can be shared with the patient directly. Alternatively,

the workstation can be coupled to a color printer. The user would simply print out hard

copies of the screen shots showing the virtual patient model.

The manner in which the software is written to provide tools allowing for simulation

of various parameters can vary widely and is not considered especially critical. One

possibility is a Windows-based system in which a series of icons are displayed, each icon

associated with a parameter. The user clicks on the icon, and a set of windows are displayed

allowing the user to enter new information directing a change in some aspect of the model.

The tools could also include slide bars, or other features that are displayed to the user and tied

to specific features of the patient's anatomy. Treatment planning icons for moving teeth are

disclosed in the published PCT application of OraMetrix, Inc., W0 01/80761, which gives

36

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

10

15

some idea of the types of icons and graphical user interface tools that could be used directly

or adapted to simulate various parameters.

Once the user has modified the virtual patient model to achieve the patient's desired

feature and smile, it is possible to automatically back-solve for the teeth, jaw and skull

movement or correction necessary to achieve this result. In particular, the tooth movement

necessary can be determined by isolating the teeth in the virtual patient model, treating this

tooth finish position as the final position in the interactive treatment planning described in the

published OraMetrix PCT application, W0 01/80761, designing the bracket placement and

virtual arch wire necessary to move teeth to that position, and then fabricating the wire and

bracket placement trays, templates or jigs to correctly place the brackets at the desired

location. The desired jaw movement can be determined by comparing the jaw position in

the virtual patient model's finish position with the jaw position in the virtual patient model in

the original condition, and using various implant devices or surgical techniques to change the

shape or position of the jaw to achieve the desired position.

The virtual patient model as described herein provides a common set of data that is

useable in a variety of orthodontic or other treatment regimes. For example, the initial and

final (target) digital data sets of the patient's tooth positions can be relayed to a manufacturer

of customized transparent removable aligning shells for manufacture of a series of aligning

devices, as taught in the Chisti et al. patents cited previously. Alternatively, the tooth

positions may be used to derive customized bracket prescriptions for use with a flat planar

archwire.

10

15

20

The choice of which treatment modality, and whether to use any additional treatment

or therapeutic approaches (including surgery) will depend on the patient in consultation with

The integrated environment proposed herein provides essentially a the treating physician.

37

platform for a variety of possible treatment regimes. Further, the creation and display of the

virtual patient model provides for new opportunities in patient diagnosis and sharing of

patient information across multiple specialties in real time over communications networks.

Figure 19 is a block diagram of an integrated workstation environment for creation of

the virtual patient model and diagnosis, treatment planning and delivery of therapeutics.

The workstation environment shown in block diagram form in Figure 19 may incorporate

many of the hardware aspects shown in Figure 1, including scanning or imaging devices

28/36 for capturing two dimensional images, such as a color digital camera or X-Ray

machine. The workstation environment will preferably include scanning or imaging devices

30/36 for capturing three dimensional images and creating 3D models of the patient,

including one or more of the following: laser scanners for scanning a plaster model of the

teeth, optical scanner such as the OraMetrix hand-held scanner 30 referenced in Figure 1, CT

scanner or MRI. In some instances, the scanning devices may be located at other facilities,

in which case the 3D scans are obtained at another location and the 3D data is supplied to the

workstation 10 (Figure 1) over a suitable communications channel (Internet) or via a disk sent

in the mail.

5

10

15

20

The workstation includes a memory storing machine readable instructions comprising

an integrated treatment planning and model manipulation software program indicated

generally at 300. The treatment planning instructions 300 will be described in further detail

below. The treatment planning software uses additional software modules. A patient history

module 302 contains user interface screens and appropriate prompts to obtain and record a

complete patient medical and dental history, along with pertinent demographic data for the

38

patient.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

A module 304 contains instructions for designing custom dental and orthodontic

appliances. These appliances could include both fixed appliances, e.g., brackets, bands,

archwires, crowns and bridges, surgical splints, surgical archwires, surgical fixation plates,

laminates, implants, as well as removable appliances including aligning shells, retainers and

partial or full dentures. In one possible embodiment, the module 304 may be located and

executed at the site of a vendor of custom orthodontic applicants. The vendor would receive

an order for a custom appliance specifically to fit an individual patient. Module 34 would

process this order and containing instruction for designing the appliance to fit the individual

morphology and condition of the patient. The vendor would take the appliance design,

manufacture the appliance in accordance with the design, and then ship the custom appliance

to the practitioner. Examples of how the appliance design module 304 might be

implemented include the appliance design software developed by OraMetrix and described in

the published PCT patent application cited previously, the customized bracket, jig and wire

appliance design software of Ormco described in the issued Andreiko patents U.S. Patent

(see, e.g., 5,431,562) and in the published patent application of Chapoulaud, US patent

publication no. 2002/002841, the techniques of Chisti et al., U.S. Patents 6,227,850 and

6,217,325, all incorporated by reference herein.

The treatment planning software 300 also obtains information on standard ("off the

shelf") dental or appliances from a module 306, which stores manufacturer catalogs of such

appliances, including 3D virtual models of the individual appliances.

The treatment planning software includes a module 308 that allows the user to input

selections as to variable parameters that affect the visual appearance of the patient, as input to

a craniofacial analysis module 328 described below. The variable parameters include patient

factors: age, weight, sex, facial attributes (smile, frown, etc.). The variable parameters also

39

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

5

10

15

include parameters affecting the teeth, including texture (color), position, spacing, occlusion,

The variable parameters further include various illumination parameters, including

reflectivity of the skin, ambient light intensity, and light direction. These parameters are

accessed though appropriate icons on the screen display, such as the icons shown in Figures

4-7, and pop-up displays that appear that prompt the user to enter or vary the selected

variable parameter.

10

15

20

The treatment planning software further uses a diagnosis and simulation module 310

that displays diagnosis data graphically and/or in report format. This diagnosis data includes

teeth position, 3D face and smile appearance, and various facial attributes.

The software further includes third party practice management software 312.

Information about treatment planes generated by the craniofacial analysis module 328 is

input to the practice management software 312. Based on the treatment plan, this software

generates the most productive scheduling of appointments for the patient. The practice

management software 312 also generates reports, provides insurance and benefit tracking,

and supports electronic claims filing with the patient's insurance company. Preferably, the

practice management software provides a flexible scheduling of patient appointments based

on progress of treatment of the patient's craniofacial anatomical structures. The progress of

treatment is obtained from periodically obtaining updated three-dimensional information

regarding the progress of treatment of the craniofacial features of the patient. For example,

the patient is periodically rescanned during the course of treatment. A new virtual patient

model is created. Depending on the progress of treatment (e.g., movement of the teeth to

target positions) the patient may be scheduled for more or less frequent visits depending on

40

their progress.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive

(312) 913-0001

Chicago, IL 60606

Referring again generally to the treatment planning software 300, the software

includes a 3D model generation module 314 that uses as input the 2D and 3D scanning

devices. A 3D virtual model of the patient is created by module 314, for example, in the

manner described previously in Figures 2 and 3.

The system further includes a custom appliance management module 315. This

module provides appliance specifications and 3D geometry data to the vendor site for the

purpose of providing necessary input for the design and manufacture of custom appliances,

such as custom orthodontic appliances, for the patient. This module also provides updates to

an appliance data module 324 for storing custom appliance data within the database. The

module 324 is responsible for managing the database of all the appliances, including custom

appliances.

5

10

15

20

The 3D virtual patient model is supplied to a knowledge database 316.

knowledge database includes a 3D Geometry data file 316 that stores the 3D geometry data

of the virtual patient model. This data is supplied to a tagging module 322 and a morphable

model module 320. The morphable model module 320 includes instructions for creating a

morphable model from various 3D model samples, using the techniques for example set forth

in the article of Blantz et al., A Morphable Model for The Synthesis of 3D Faces, Computer

Graphics Proceedings SIGGRAPH '99 (August, 1999).

The tagging module 322 includes instructions for tagging or placing pieces of

information regarding the virtual patient model into each patient record, which is used for

statistical procedures. In particular, the tagged information is supplied to a meta-analysis

module 326. The meta-analysis module implements a set of statistical procedures designed to

accumulate experimental and correlational results across independent studies that address a

related set of research questions. Meta-analysis uses the summary statistics from individual

41

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

studies as the data points. A key assumption of this analysis is that each study provides a

different estimate of the underlying relationship. By accumulating results across studies, one

can gain a more accurate representation of the relation than is provided by the individual

study estimators. In one example, the software will use previous patient cases/studies to help

in the craniofacial analysis module 328. For example, surgery cases for "lip and chin" will be

one set of independent studies, whereas jaw surgery to correctly position the upper and lower

jaw will be another. An orthodontist trying to align the upper and lower jaw will do a meta-

analysis with the module 326 in order to see how this treatment will affect the patient's lip

and chin.

10

15

20

The output of the morphable model from module 320 and the meta-analysis from

module 326 is provided to a craniofacial analysis module 328. This module takes as input,

patient information and the patient 3D virtual model to generate diagnosis and simulation

data. Based on one or more simulation results, this module 328, and/or module 330

generates a treatment plan and appliance selection. User involvement is contemplated in

modules 328 and 330. In particular, the user may interact with the patient information and

the morphable model, and vary the parameters 308, to simulate different possible treatments

and outcomes to arrive at a final or target treatment objective for the patient. The

craniofacial analysis module 328 may include some or all of the treatment planning features

described at length in the published PCT application of OraMetrix, Inc. cited previously.

The software instructions included in the craniofacial analysis module 326 preferably

includes a set of instructions providing the user with user interface tools (e.g., icons), for

visually studying on the user interface 316 the interaction of the craniofacial anatomical

structures and their relationship to the external, visual appearance of the patient. For

example, tools may provide a chewing simulation. Alternatively, the tools may provide a

42

smile function in which the face is morphed to smile, showing the position of the teeth, gums,

lips and other structures. These tools simulate changes in the anatomical position or shape of

craniofacial anatomical structures (teeth, lips, skin, etc.) and show the effect of such changes

on the visual appearance of the patient. As another example, the tools may include tools for

modifying the shape or position of one or more bones of the upper and lower jaws, and show

how those modifications affect the patient's appearance and smile.

With reference to Figure 7, the user would activate one of the icons 35 at the top of

the screen display. The icon may be associated with a function that would allow the user to

reposition the location of the upper and lower teeth. After the user changes the position of

the teeth, the user would activate another icon, "smile", and the face would morph to a smile

with the teeth in the new position.

After the patient simulations have been completed and the patient and physician are

satisfied, the resulting data set of teeth position, jaw position, etc. are stored by the diagnosis

and simulation module 310 of Figure 19. This module 310 preferably includes a routine for

storing a three-dimensional representation of said patient's craniofacial structures (e.g., teeth)

in a format suitable for use by a manufacturer of orthodontic appliances. Each manufacturer

may have a unique format needed for use by the manufacturer, and the routine takes that into

consideration in storing the data. For example, a manufacturer may require 3D digital models

of the teeth in initial and final positions in the form of triangle surfaces, along with archwire

and bracket prescription data.

It is contemplated that the creation and usage of the virtual model may occur at the

patient care site. In particular, the treating physician or orthodontist will access the scan and

photographic data, create the virtual model therefrom, and perform the treatment planning

and simulation described herein in their own office. Once the treatment plan is arrived at,

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

15

20

the treating physician can export the virtual patient model or some subset of data to appliance

manufacturers or specialists, as indicated in Figure 1.

Alternatively, the virtual patient model may be created at a remote location. In this

latter example, a third party, such as an appliance manufacturer, may be the entity that creates

the virtual patient model and makes it available to the treating physician. In this example, the

treating physician will have access to the scanners, X-Ray, digital camera, or other imaging

device, obtain the required data from the patient, and forward such data to the third party.

The third party executes the instructions to create, visualize and manipulate the virtual patient

model. This model can be transmitted to the treating physician for their review and usage.

Then, either the third party could create a proposed treatment for review and approval by the

treating physician, or the treating physician could create the treatment plan. The plan is then

transmitted to one or more appliance manufacturers for fabrication of therapeutic devices

(e.g., brackets and wires, aligning shells, maxillary expansion devices, etc.)

A treatment plan created from the virtual patient model described herein may be one

in which only one type of appliances, e.g. fixed of removable, is used during the entire course

of the treatment. For example, the treatment plan may be one in which brackets and wires

are the type of appliance that is used. Or, alternatively, the treatment plan may be one in

which removable aligning shells are the type of appliance that is used.

On the other hand, the treatment plan might be such that it is a hybrid plan requiring

the use of different types of appliances during the course of the treatment. In the hybrid

orthodontic treatment plan, a variety of scenarios are possible. In one type of hybrid

treatment plan, different types of appliances might be used at different times during the

course of the treatment. For example, patient may start out with brackets and wires and shift

at some point during treatment to an approach based on removable aligning shells.

44

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

10

15

another type of hybrid treatment plan, different types of appliances might be used

simultaneously, for example in different portions of the mouth, for example brackets and

wires could be used for certain teeth and transparent aligning shells uses for a different set of

teeth. A hybrid treatment plan may be chosen right from the beginning, or it may be

introduced dynamically at any stage during the treatment course.

To develop a hybrid treatment plan, the treatment planning software will preferably

include features of the appliance design and treatment planning software of the manufacturers

of the appliances that are used in the hybrid treatment. As one example, the treatment

planning software may include the wire and bracket features of the OraMetrix treatment

planning software described in the published application WO 01/80761, as well as the

treatment planning software described in the Align Technologies patents to Chisti et al., U.S.

Patents 5,975,893 and 6,227,850. The software would thus allow the user to simulate

treatment with brackets and wires for part of the tooth movement to reach a particular

milestone, and also design the configuration of intermediate tooth positions and configuration

of removable aligning shells for the remainder of tooth movement. Alternatively, the shape

of the aligning shells could be determined automatically via the treatment planning software

from the tooth configuration at which the shells are first introduced to the patient and the final

tooth position in accordance with the teachings of the Chisti et al. patents.

Figure 20 is an illustration of the integration of the patient data acquisition, treatment

planning and appliance design functions that are facilitated by a preferred embodiment of the

unified workstation 14. The workstation is provided with a plurality of image data sets 400,

which can include 2D data (e.g., photographs) 402, 3D image data 404 from various 3D

image sources, static models 406 of all or part of the patient's craniofacial anatomy, dynamic

models 408 of all or part of the patient's craniofacial anatomy, color models 410, and

45

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

5

10

15

possibly other types of image data. The workstation 14 includes software 314 (such as described above in conjunction with Figure 19) that takes any possible combination of this image data to produce a virtual patient model 34. From this virtual patient model, the workstation in one possible embodiment includes one or more treatment planning tools or software 300 for planning treatment for the patient. These treatment planning tools could include specific software provided by vendors of treatment planning software or appliances, such as manufacturer #1 software 412, manufacturer #2 software 414, software for manufacturers # 3, 4, 5, . . . at 416, 418, 420, as shown. Such software would be operable on the virtual patient model 34 and associated data sets representing the teeth as described at length herein. To provide interoperability of the software on the virtual patient model, the virtual patient model may have to have representations of the data that is compatible with the software of various vendors, which is within the ability of persons skilled in this art. Moreover, once appliance designs have been created by the various species of treatment planning software, the preferred embodiment of the workstation allows export of appliance design, tooth position data or other required outputs to any appliance manufacturer so as to allow the manufacture of a customized orthodontic appliance. In other words, if the workstation is equipped with OraMetrix treatment planning software, such software could output tooth position data, appliance design data and any other required data into a format compatible with the manufacturing requirements of any appliance manufacture. interoperability of data formats for appliance design is shown by arrows 421. Thus, the workstation provides a conversion or formatting of appliance design data into a data set or output format specified by any one of a variety of particular appliance manufacturers. In the illustrated embodiment, the available therapeutics data sets are shown as manufacturer no.1 data set 422 (brackets and customized wires), manufacturer no. 2 data set 426 (brackets and

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

10

15

wires), manufacturer no. 3 data set 426 (removable aligning shells), manufacturer no. 4 data

set 428 (brackets and wires), or still other sets 430. The appliance design data set is then

furnished over the Internet to the vendor of such appliances for manufacture of a custom

appliance. Hybrid treatment plans, as described above, are one possibility of a treatment

plan that may be developed using the workstation and virtual patient model described herein.

In one possible variant of the invention, the treatment planning software tools 300 are

also provided at a remote location and some of the tasks of appliance design may be

performed as a service by a separate workstation, such as a workstation of an appliance

manufacturer. In this situation, the virtual patient model 34 could be provided to the

appliance manufacturer, a proposed treatment plan is prepared and furnished to the

practitioner, and after the plan is approved, the appliance manufacturer coordinates the

furnishing of appliance design data to any designated appliance manufacturers that are used

to furnish the custom appliance.

5

10

15

20

In one possible embodiment, the treatment planning software includes a set of

instructions that perform a measurement function to measure distances in two or three

dimensions in the virtual patient model, e.g., arch form shape measurements, and compare the

measurements with reference dimensions for an "average" patient of similar age, sex, and

race. These measurements could be obtained in any convenient manner, for example from

textbooks, organizations, practitioners, etc. These measurement tools would be invoked

during the course of treatment to compare tooth movement and current tooth position with

expected positions and if deviations occur, the variances could be used as information to

47

modify one or more aspects of the treatment plan, such as change the appliance design.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

A specific example of interactive treatment planning with a unified workstation will

be now described in conjunction with Figures 21-58. Figures 60-91 describe evaluation of a

proposed set-up that has been generated in the manner described in Figures 21-58. In order

to provide the functions and features described in the features, the workstation includes a

computing platform (general purpose computer) having a graphical user interface, a processor

and a computer storage medium containing digitized records pertaining to a patient. The

digitized records include image data, such as photographs, x-rays, and scan data of the teeth.

The workstation further includes a set of software instructions providing graphical user

interface tools for providing access to the digitized records, such as display and manipulation

of the images or scan data in the form of 3D models. The workstation further includes

software instructions for execution by the processor for facilitating treatment planning for a

patient.

10

15

20

While there are various ways in which practitioner may go about the process of

designing a treatment plan for a patient, in a preferred embodiment of the invention certain

specified tools are provided which allow a treatment plan to be developed in which

constraints can be identified and the treatment plan can be developed without violation of

such constraints.

Referring now to Figure 21, a screen shot from the graphical user interface of the

workstation of Figure 1 is shown. The workstation includes a computer memory that stores,

and makes available to the practitioner, records in the form of digital data pertaining to some

or all of the following: the patient's clinical history, medical history, dental history, and

orthodontic history as well as 2D photographs, 2D radio graphic images, CT scans, 2D and

3D scanned images, ultrasonic scanned images, and in general, non-invasive and sometimes

invasive images, plus video, audio, and a variety of communication records, such notes,

48

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive

Chicago, IL 60606 (312) 913-0001

records of office visits, patient letters or communications, etc. All records and images are

digitized. The records and images are made available through suitable user interface icons

which cause display of the images and records on the user interface. The images can be

combined or superimposed to create a virtual patient model that includes surface features

(soft tissue) of the patient in one possible embodiment.

The workstation also further maintains a comprehensive set of computer instructions

providing tools in the form of icons, screen displays, windows, functions and features,

accessible through the user interface of the workstation to assist the practitioner in planning

the treatment. Various types of tools are contemplated; numerous examples are set forth

herein.

10

15

20

In Figure 21, a set of tabs 450, 452, 454, 456 and 458 are provided. The tab 450 is a

patient information tab which provides suitable screen displays for entering a variety of

patient information, such as their name and address, dental and clinical history, insurance

information, diagnostic information, names of other treating or consulting practitioners, etc.

Tab 450 will be discussed in further detail below in conjunction with Figures 21A and 21B

and Figures 54-58. The tab 452 is a tab whereby the user accesses scan or other image data

and accesses instructions and menus for scanning the patient with an in-vivo intra-oral

scanner such as described in the previously cited OraMetrix PCT application.

tab by which the user accesses the digital 3D impression of the teeth, obtained from the

scanning of the patient. Tab 456 is a case management tab and includes a number of specific

available screen displays and menus which are shown in the menu 460 in the lower left of the

Figure. The case management tab, and its various features, is described at length in the

following discussion. Additionally, there is a digital treatment planning tab 458 which

provides further menus, tools and displays by which the practitioner may further move teeth

49

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

and design the shape and configuration of a customized orthodontic appliance. An example

of the types of menus and tools that are available in the tab 458 is the OraMetrix treatment

planning software described in application serial no. 09/834,412, filed April 13, 2001.

However, it is possible to provide, in the workstation, a suite of treatment planning software

from different appliance manufacturers in which case the user could access the treatment

planning software for whatever appliance manufacturer the practitioner wished to use for

treatment of the patient. In this situation, it may be necessary to format the tooth data in a

format compatible with the appliance design and treatment planning software so as to ensure

compatibility between the various systems that may be installed on the workstation.

In Figure 21, the user has selected a "treatment strategy" icon 461, which causes the

display 462 to appear. In this display, there is a field 464 for the user to enter high level

diagnosis and problem classification information, for example in the form of text. A field

466 is provided which provides a matrix format by which the conditions relevant to the

patient's soft tissue, skeletal, and dental anatomy are entered, each with respect to vertical,

sagittal, and transverse positions, again in text form. The display also includes a treatment

strategy field 468 where the user will indicate the general, high level approach to treatment,

such as any proposed extractions, appliance type, stages of treatment, etc. These fields 464,

466 and 468, along with displayed image data for the patient, assist the practitioner in

identifying the constraints pertinent to the treatment planning.

Figure 21A shows the patient information tab 450, with the slide bar 850 moved next

to "history". The screen display shown in Figure 21A appears, with a field 451 for which

the user can enter alerting medical conditions, such as AIDS or HIV infection, epilepsy,

allergy conditions, tuberculosis, etc., along with any explanatory data or comment. In field

453, the user is provided with tools to enter general medical condition regarding the patient

50

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

10

15

by clicking on the drop-down menu as shown, and entering any appropriate data or

commentary, as shown.

In Figure 21B, the user has moved the slide bar 850 to the "Examination" icon, which

causes the display shown in Figure 21B to appear. This screen allows the user to enter

dental examination data in field 455, a tooth chart 457 where the user clicks on a particular

tooth and enters tooth data in fields 459, 461 and 463 as indicated.

After the user has entered the information into the fields 464, 466, 488 shown in

Figure 21, the user clicks on one of the other icons in the field 460 to continue the process of

case management and initial treatment planning. At this point the information entered into

the fields of Figure 21 is stored in the computer memory of the workstation.

In Figure 22, the user has now selected a "Midline and Aesthetic Occlusal Plane" icon

470, which causes the screen display 472 to appear. The user uses this screen to evaluate

and define both vertical and horizontal lines of references such as soft tissue midline,

interpupilliary line, etc and also define the dental midlines for the upper and lower dentition

and the aesthetic occlusal planes for both the upper and lower arches and cant of the occlusal

planes These midlines and occlusal planes are designed relative to the face, which here is a

global reference. These lines are useful references for where you want the patient's teeth to

move.

10

15

20

When screen display 472 is activated, the workstation displays a pair of two-

dimensional color photographs of the patient, shown as a photo 474 with the patient's mouth

closed, and a photo 476 with the patient smiling. The display includes a field 478 where the

patient can maintain the midline that the user marks on the images, as described below, or

activate one of the other tabs indicating treat to upper midline, treat to lower midline, or

provide a custom midline. The midline is entered using the tools 486 on the right hand side

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive

Chicago, IL 60606 (312) 913-0001

of the screen A region 480 is provided for the Aesthetic Occlusal Plane (occlusal plane for

the front teeth), which the user can indicate or mark on the images of the patient using the

tools 486 on the right hand side of the screen. The user marks an Aesthetic Occlusal Plane

(AOP) for both the maxilla and mandible dentition, and the user is provided with fields 480

and 482 for customization of these planes (technically, lines in two dimensions). A tab 484

is provided to create a customized canted AOP with various tabs as shown. Thus, the tools

provide the user to mark, among other things, a midline and maxilla and mandible levels and

cant of an aesthetic occlusal plane.

10

15

20

The display of Figure 22 includes a set of tools 486 in the form of icons which, when

selected, allow the user to mark on the images various vertical and horizontal lines. For

example, the user can mark an upper occlusal plane on the photographs of the upper arch of

the patient, a lower occlusal plane (line) in the lower arch of the patient, and marking

various positions in the upper and lower occlusal planes, e.g., marking a posterior position of

the upper or lower occlusal plane (line in 2D); marking a functional position of the upper or

lower occlusal plane; and marking an aesthetic position of the upper or lower occlusal plane.

As shown in Figure 23, when the user activates the legend "L" icon 488, a window

490 pops up and a legend appears that explains the acronyms for various lines and midlines

that the user may mark on the images.

As shown in Figure 24, the user has activated various icons 486 and has drawn on the

virtual model of the patient an aesthetic upper occlusal plane ("AU") 494 and a aesthetic

upper perpendicular line ("AUP") 492 in the left-hand image, and a treatment upper occlusal

plane ("TxU") 498 and a treatment upper perpendicular line ("TxUP") 496. The lines 492,

494, 496 and 498 are all user specified in terms of their location. The location is selected by

using the workstation mouse, moving the cursor to the location where the user wishes to draw

52

the midlines and occlusal planes, and clicking the mouse.

Referring now to Figure 25, after the user has proceeded to mark the location of midlines and occlusal planes, the user clicks on the icon 500, at which point the "Occlusal Plane and AP [Anterior, Posterior] Positions" screen display 502 appears. The display 504 shows other aspects of the patient virtual model, here a two-dimensional lateral x-ray view of the face, including teeth and jaws. Shown on the display 504 are two of the occlusal planes (line in 2 dimensions) the user has previously entered from the screen display of Figure 24: a normal occlusal plane ("NO") 506, and a Treatment Occlusal Plane Upper (TxOU) 508, indicating the occlusal plane to which the teeth are designed to be aligned with as a result of treatment of the patient. These lines can be segmented into three separate lines, one for the posterior, functional and aesthetic. The screen display includes a set of icons 510 providing tools for the user to mark various occlusal planes and locations thereof in two dimensions, as well as a Legend icon 520. The display also includes a region 512 whereby the user can modify the location of the aesthetic occlusal plane (front teeth, teeth 1-3), a functional occlusal plane (teeth 3-6), and a posterior occlusal plane (rear teeth, teeth 6-8), for both the upper and lower jaws. The axis of cant of the occlusal plane can be changed by rotating around a predetermined center or fulcrum of rotation. Also, the A/P position and inclinations and the vertical relations of the incisors with respect to the occlusal plane can be represented by animating teeth as shown in Figure 56, 57. The desired position of the incisors can be planned. The position of the incisors also drives the change in the position of the soft tissue, e.g. lips. Any changes can be compared against a normative database of positions of various craniofacial structures. A complete 2D cephalometric analysis is thus possible.

As shown in Figure 26, when the activates the legend icon 520, the window 522 pops up and provides a legend for the acronyms which accompany the occlusal planes that are

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

10

15

shown on the x-ray image 504. The various occlusal planes 522 are accessed by activating

the icons 510 at the right hand side of the screen display, and using the mouse to indicate

their location on the image 504.

5

10

15

20

Referring to Figure 27, the user has finished marking the midline, the occlusal

plane(s) and making any anterior/posterior adjustments, and proceeded activate the

"Reference Tooth" icon 550. This action causes the display 552 to appear. The display 552

includes a reference tooth selection field 556. The display also changes to show a panoramic

X-ray 554 to appear, showing all of the patient's teeth to appear in a line, along with

associated soft tissue and surrounding bone structures. The user makes a judgment decision

as to which tooth or teeth should be moved the least (or not at all), and selects this tooth (or

teeth) to be the reference tooth or teeth. This action is completed by the user moving the

mouse cursor to the tooth or teeth in the field 556 and clicking the tooth they wish to select as

the reference tooth. Here, the user has selected tooth 3 on the upper right hand side as the

reference tooth, as indicated at 558. Any changes in crown position in two dimensions or

root positions are seen and transferred into the three-dimensional data sets.

Figure 27A shows a screen display with the display showing a flattened 3D x-ray of

the teeth flattened such that all the teeth lie in a two dimensional line, with each tooth having

a tooth axis indicated by a line segment having one end terminating at the cusp or tip of the

tooth and the other end of the line segment terminating at the end of the tooth root. This

screen is displayed simultaneously with the field 556 showing the reference tooth selected.

The user, having inspected the axes of the teeth and their relationship to other teeth via the X-

ray, may select a different reference tooth by simply clicking on a different tooth in the field

556. Typically, the user will select a reference tooth in which the axis of the tooth does not

move during the course of treatment, and the displays of Figure 27 and Figure 27A facilitates

54

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

that selection. The screen display of Figure 27A facilitates the measurement of the tooth axes

to make the reference tooth selection.

The treatment planning process continues by using the graphical user interface to

align two dimensional images of the patient, e.g. x-rays, with three-dimensional virtual teeth

models. In this manner, the user progresses from two-dimensional treatment planning to

three-dimensional treatment planning. One possible embodiment is shown in Figures 28-31.

In Figure 28, the user has selected the "Alignment 2D-3D" icon 570, which causes the screen

display 572 to appear. In this display, a 3D virtual model of the teeth 116 appears on the

display, superimposed over a two dimensional X-ray photograph 505. The 3D model of teeth

116 is created by suitable scanning techniques, such as in-vivo scan of the dentition or a scan

of a model of the dentition, as described previously. After the scan is obtained, the teeth are

separated from surrounding structures and represented in the computer as individual,

individually moveable, virtual tooth objects. The display includes navigation and other

icons by which the user can rotate the model 116 in any desired orientation, show only one or

both arches, select or deselect for display gingival tissue, occlusal planes, etc. The 3D tooth

models are scaled so that they coincide in size with the size of the teeth in the 2D image. The

superposition shown in Figure 28 could be either performed manually, or possibly

automatically with suitable pattern matching algorithms to identify tooth objects in the 2D

image and align the 3D tooth objects with the teeth in the 2D image.

In Figure 28, the functional occlusal plane 508 is displayed together with the teeth and

the x-ray. Whereas in Figure 24, the upper occlusal plane was shown as merely a line, in

Figure 28 the occlusal plane 508 is represented in two dimensions but it actually is also

represented in three dimensions in Figure 43. Thus, the original 2D representation is

transferred to a surface in three dimensions. The user is able to view the arragement of the

55

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

10

15

teeth relative to the bone in any orientation, such as front perspective or side perspective.

The arrangement in Figure 28 facilitates the user understanding the relationship of the

3D teeth with respect to the soft tissues as well as bone. The 3D plan of the teeth can be

oriented relative to the occlusal plane 580 that the user has defined.

In Figure 28A, the user has activated icons 574 for display of both arches and a

midline plane 581.

5

10

15

20

In Figure 29, the user has selected for display both the upper and lower arches by

activation of icon 602, and treatment occlusal planes 606 and 608 for the upper and lower

arches. The occlusal planes 606 and 608 are activated by selecting the icon 604. Note that

in Figure 29, the occlusal planes are rendered and displayed in three dimensions as a three-

dimensional surface, whereas, initially in Figure 24, the midline and occlusal planes were

rendered in two dimensions. The three dimensionality of the planes 606 and 608 is hard to

see in Figure 29, but becomes more apparent when the model of the teeth is rotated or viewed

from an orientation that is not so closely in line with the planes 606 and 608.

In Figure 30, the user has activated icon 610 which causes a mid-sagittal clipping

plane 612 to appear. The location where the clipping plane 612 intersects the front teeth in

the upper and lower arches is shown by the shaded areas 614. The clipping plane can be

moved over the arch to view the teeth in any cross-sectional view, using navigation icons.

The view shown in Figure 30 allows the user to judge subjectively the relationship between

the upper and lower incisors, and compare that with the 2D views, for example, from a 2D X-

ray.

As shown in Figure 31, the user can adjust the position of the clipping plane and

thereby change the location at which the plane intersects the upper and lower arches, as will

be appreciated from a comparison of the shaded areas 614 in Figure 31 with the areas 614 in

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

Figure 30. The user adjusts the position of the upper and lower incisors in the clipping plane

illustration to match the position that was determined in the 2D lateral view in Figure 57.

When the user is satisfied with the 2D - 3D aligning step, the user proceeds to

additional tasks, including space evaluation, and space management tasks by which the user

first evaluates how much space is needed to align the teeth and then how he wishes to

manage the space. The user further proceeds with the design of an desired arch form.

This is done for both arches, typically the mandible first and then the maxilla. However, at

any time the user can view both arches together by activating a hide/display icon. To proceed

to these tasks, in the illustrated embodiment, the user selects a mandible space management

icon 620, as shown in Figure 32, which causes the screen display 622 to appear. The screen

display 622 includes a plurality of icons 624 which are used for hiding and displaying

various aspects of the virtual patient model, soft tissue, occlusal planes, and other features of

The central region of the display 622 is used to display the 3D virtual teeth the software.

116. The display 622 also includes a lower region 626, where the user can activate an Arch

Form tab 636, as shown, or other tabs, including a Vertical Positions tab 628, a A/P positions

tab 630, and an Extractions and Stripping tab 632. The arch form tab 634 includes an area

indicating that the user has selected a customized arch form. However, by activating the drop

down icon 638, the user can scroll through and select pre-defined arch form types that are

stored in the workstation, and adapt tooth position to standard arch forms. At all times, the

user is able to interactively move any tooth or teeth on the graphical user interface relative to

the desired arch form, by clicking the tooth to select the tooth and then dragging the tooth

with the mouse to a new position.

Space analysis can be dynamically evaluated by affecting the following parameters:

midline, arch form, A/P position, tooth position, the reference tooth, tooth size, spatial

57

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

5

10

15

distribution of the teeth in the arch and by appliance prescription, either selectively or in

tandem. Furthermore, space management can be effectuated by simulation of interproximal

reduction, buildup of the tooth, extraction, distal and mesial tooth movement, expansion of

the jaw, axial inclination angle change, rotation change, overjet and overbite change,

appliance choice, adjustment of inter-arch relationship, or selectively maintaining crowding.

The tab 634 further includes measurement tools 640 which provide cuspid distance

measurements and inter-molar distance measurements for the current tooth positions

displayed on the screen. The user can also set points anywhere on the virtual model and

activate an icon to get a distance measurement, or invoke a graph tool as described elsewhere.

Figure 32 also shows a display 644 that provides an arch length discrepancy measurement (in

terms of mm) which indicates, given the current virtual tooth positions, whether there is

sufficient length in the arch (positive values) or whether some interproximal reduction, tooth

rotation, extraction, distal movement of the molars, uprighting of the molars, changing of the

torque of the teeth, changing the A/P position of the incisors, expanding the arch form,

maintaining selective crowding, adjusting the level of the occlusal plane or the midline, axial

inclination of teeth, overjet or overbite or other action is required to fit the teeth in the arch

(negative values). The left-right balance in the arch length discrepancy can be changed by

interactively selecting the midline and moving the midline either to the right or left.

In Figure 32, the presenting malocclusion is seen. The user has indicated that they

wish to hold the malocclusion as indicated at 646. In Figure 33, the keep lateral segments

tab 650 has been activated. This feature immobilizes the lateral segments (premolars and

molars) and allows for the alignment to occur only through the change in the arch form

anteriorly. The user can define the incisor position based on the lower incisor based on

where it was initially. The user can hold the buccal segment (basically holding the lateral

58

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

10

15

segments) then impose no movement and you allow for full alignment of the anterior teeth.

What this has done is tell the practitioner how much space is needed in the arch length and at

the same time the display provides the inter-molar and inter-cuspid width based upon holding

In other words, as the user changes the constraints in terms of the posterior teeth fixed.

tooth mobility, the user is provided instantaneous measurements in terms of arch length and

the new positions of the teeth that are allowed to move. The user can selectively immobilize

teeth or allow their free movement in three planes of space, either individually or in groups.

The teeth are moved to a more ideal position. This action changed the cuspid distance

from Figure 32, but did not change the molar distance. The length discrepancy tool 644

indicates that the right and left arch length has increased from the previous values, which

would require some interproximal reduction, extraction, or some other action to fit the teeth

to the arch length.

10

15

20

In Figure 34, the user has selected "full aligned", as indicated at 660.

cuspid distance and the molar distance have changed, as indicated by the measurement tools

The length discrepancy, as indicated by the display 644, has now returned to zero,

indicating that the arch form is "good". "Fully aligned" allows free movement of all the teeth

along the pre-selected arch form. If some arch length discrepancy remained, the user could

simulate a change in arch form, midline, occlusal plane, rotation of teeth about their axis,

extraction etc. to rectify the situation. If the user wants to customize the shape of the arch

form they activate the slide line tab, discussed later.

In Figure 35, the user has activated icons 670 and 672, which causes the display to

show both the original position of the teeth (areas with dark shading indicated at 674) and the

new position as a result of the mandible space management exercise (white areas, indicated at

676. This color coding helps the user visualize the tooth movement that will occur in

59

accordance with the proposed treatment plan. This feature of reference back to the original

malocclusion is available at any time and in any plane, both in 3D or 2D images or

combination.

10

15

20

In Figure 36, the user has activated icon 680 which causes a slide line 682 to appear.

The slide line can positioned at different levels, such as the bracket level, the level of the cusp

tips, or the level of the interproximal contact points.

Arch length discrepancy can be defined at various levels, including contact points,

cusp tips and at the level of brackets, based upon the location of the slide line that is chosen.

Then, the effect of bracket prescription on the dentition is also modeled in defining the

position of the teeth in the arch, thus providing the clinician with a method of understanding

the effects of his appliances on arch length inadequacy.

The slide line 682 is a tool that assists the user in changing the shape of the arch form.

The slide line 680 includes anchor points 683 spaced along the length of the slide line 682,

which are affixed to labial surfaces of the teeth in the positions shown, The slide line 682 also

includes points 681 equidistantly spaced from the anchor points, which the user manipulates

to cause the slide line to bow out or in relative to the teeth, and thereby change the shape of

the arch form. For example the user would click on one of the points 681 and drag the point

681 out away from the slide line, which would cause the slide line to bow outwardly towards

the point 681. The clamping or anchor points can be moved by the user anywhere along the

slide line. The slide line (as was the case with the midline) allows for designing asymmetric

arch forms. Whenever the user wishes to compare the proposed arch form with the original

tooth position, they activate an icon at the top of the screen and the original tooth position is

60

also shown, with the difference in position shown in a contrasting color.

Figure 36A shows a screen display when the treatment planning icon 900 is actuated.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

This icon takes the user to addition treatment planning software which enables a user to

further define the shape of the arch form, move teeth relative to each other, and design

customized appliances. Many of the features shown in the display of Figure 36 A are

described in the published PCT application WO 01/8076, and therefore will not be repeated

here. The user has highlighted the "TARGET PARAMETER" icon 902, which allows the

user to customize the configuration of the desired target arch form. Here, in area 906, the

user can scroll through a series of available general arch form types. In field 903, the user is

provided with various tools to make the arch form symmetric (or not symmetric), enter in

values for transverse, sagittal or vertical movement, or change the relation of the overbite or

overjet. As shown in Figure 36A, the slide line feature 682 is also displayed, which the user

can use as described previously to change the shape of the arch using the slide line features.

In Figure 37, the user has activated an icon to show a midline plane 690 to appear.

By interactively moving the midline (using click and drag technique), the user can change the

orientation of the midline relative to the teeth in the arch and thereby change the right/left

values shown in the length discrepancy icon 644. The midline is shifted slightly to the

patient's right. The lateral segments are held fixed while the midline shift is only carried out

by the anterior teeth. If "full aligned" had been checked the all the teeth would move along

the slide line to accommodate the new position of the midline. Note also in this Figure that

the new position of the teeth can be measured against the original as indicated by the lightly

speckled pattern of the posterior teeth, which indicates little movement of the lateral teeth

from the malocclusion position (the difference in position indicated by the dark or shaded

spots on the teeth).

10

15

20

In Figure 37A, the user has selected a parameter tab 690, which includes icons 691 for

midline shift (left or right), the measurement tools providing cuspid and intermolar distance,

61

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

and a tool 694 indicating the amount of over bite and overjet, given the current configuration.

This data is obtainable since the position of the upper arch relative to the lower arch is known

in the in the computer, and the position of the incisors in both arches is known in three-

dimensional space. Thus, the display allows the user to interactively move the teeth and

immediately realize the effect of tooth movement on overbite and overjet. You can also

adjust overjet and overbite parameters to see their effect on the arch length inadequacy.

By activating icon 692, the user can manage the spacing between teeth by having all

spacing between teeth to be equal. By activating icon 693, the user invokes a collision

detection algorithm that prevents the user from moving teeth in a manner such that a tooth

collides with another tooth, either in the same arch or in the opposing arch. The software

allows for interproximal reduction by morphing the tooth shape to match the available space,

using a simple morphing algorithm that shrinks the tooth in two or three dimensions.

In Figure 38 the user has invoked the simulation of the midline plane 690. By

movement of the midline plane left or right they can correct the length discrepancy (left or

right) for the arch.

10

15

20

In Figures 39 and 39A, the user has activated a gingival display icon 700 which

causes the display to show both the 3D virtual teeth 116 and a 3D model of the gingival tissue

702. This icon can be activated at any time in the space planning screen. The model 702 of

the gingival tissue can be obtained from a scan of the patient's mouth with an intra-oral 3D

scanner as described previously, or it could be obtained from a scan of a physical model of

the patient. The gingival tissue 702 is separated from the scan data of the teeth 116 so that

the two types of objects can be modeled as separate 3D objects and thus displayed either

alone or together. In Figure 38, the proposed treatment indicates substantial bumps 703 on

the lingual side of the teeth. The gingival tissue is morphed to show the effect of tooth

62

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

movement to this position on the gingival tissue. Since gingival tissue closely follows the

underlying bone, the bumps 703 may indicate that the proposed tooth movement is

incompatible with the patients mandible bone structure. This could be confirmed for example

by zooming in on the location of the bumps 703, invoking the display of underlying bone

structure from an X-ray or CT scan stored in the workstation, and examining the

superposition of the virtual tooth over the underlying bone. In Figure 39A, the user has

moved tooth 701 slightly and the gingival tissue 702 "morphs" to follow the position of the

tooth.

10

15

20

In Figure 40, the user has activated navigation icons on the display to rotate the model

of teeth + gingival tissue to a new orientation. The user has also activated the A/P positions

tab 630, which allows the user to constrain either a singular movement of teeth or a group

movement of teeth in the A/P direction. The display shows a field 710 which allows the user

to either maintain or interactively change the A/P position of the right molar, a field 712 that

allows the user to interactively change the A/P position of the central incisors, and a field

714 that allows the user to interactively change or maintain the A/P position of the left molar.

As shown in Figure 40, the user has interactively changed the anterior position of the central

teeth by an amount of 1.5 mm, and the change is simulated by moving the central teeth

anterior wise by 1.5 mm and keeping the back teeth fixed and showing the new position on

the screen display. The 1.5 mm movement may not be sufficient given the 0.9 mm arch

length inadequacy remaining.

In Figure 41, the user has selected an icon 720 that causes the display to show an

occlusal plane 722, along with the position of the teeth 116 in a proposed arrangement. The

areas where the occlusal plane intersects the teeth 116 are shown as white areas 724 of the

63

teeth.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

As shown in Figure 42, the user has re-oriented the teeth to a side view and activated

the icons to show both the teeth 116 and the treatment occlusal plane 722. This view gives

the user a better view of the points of intersection of the occlusal plane with the teeth.

In Figure 43, the user has activated the extractions/stripping tab 632. This causes the

display to show the teeth of the lower arch, and a row of numbers corresponding to each tooth

in the lower arch. The display also includes a check box 730 below each number. If the

user un-checks the box 730, the corresponding tooth disappears to simulate an extraction.

The numbers above the check boxes indicate the spacing needed in each half of the arch to

meet the arch length requirements and reduce the arch length discrepancy to zero.

The display also includes an icon 721, which, when activated, all the teeth in the arch

are moved in three dimensions such that they just touch the occlusal plane. This is an

automatic function, since the location of the teeth and teeth cusps are known in three

dimensions, and the treatment occlusal plane is also known in three dimensions.

Figures 44 and 44A show the user activating a grid icon 740, which causes the display

to show a grid 742. The grid, which is in units of millimeters, gives the user an additional tool

to gauge distances between teeth (e.g., inter-molar width, inter-canine width, etc.). In Figure

44, the numbers in the field 632 above the check boxes are essentially evenly distributed

among the teeth in the left hand side of the arch, and in the right hand side of the arch,

indicating that the interproximal reduction is symmetrically and evenly distributed between

the teeth. The teeth can be morphed accordingly to simulate the change in shape either

through reduction in size or buildup in size as necessitated by the treatment. By un-checking

any of the boxes in field 632, the user can simulate a tooth extraction and move teeth to

manage the gap in space on a tooth by tooth basis, or have the space distributed evenly

between the teeth, etc. Once the user is satisfied, they can go to the screen of Figure 28 and

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

10

15

20

see where the teeth are relevant to the bone.

Figure 45 shows the virtual model of the teeth of the lower jaw. Because the teeth

are shown in a plan view it may be difficult for the user to see the texture and surface

characteristics of the teeth. Thus, the user interface preferably includes an icon 744, which

when activated allows the user to change the simulated illumination of the teeth, i.e., become

brighter, or come from a different place, in order to more clearly show the surface

characteristics of the teeth.

10

15

20

Figure 46 shows the user activating the vertical position tab 628. This tab includes a

field 750 in which the user can either maintain the malocclusion vertical position of the right

molars, or change the position by any user specified amount either towards or away from the

occlusal plane. Similarly, the tab includes a field 752 where the user can either maintain or

change the vertical position of the front teeth. A field 754 allows the user to maintain or

change the vertical position of the left molars.

Figure 47A shows a contact points feature that provides the user of a graphical

display of the contact points between teeth. The user has activated the arch form tab 634, but

the contact point function can be invoked at any time while the user is in the "Mandible

Space Management" routine. The user has activated icon 760, which signifies that the user

wishes to inspect and measure and the contact points between the teeth. When icon 760 is

clicked, points 762 appear on the teeth, which indicate the points of contact between the

teeth. The points 762 are joined by straight lines 764, which indicate the shortest distance

between the points. By placing the cursor on a particular tooth, the display shows the tooth

number and distance along the line in millimeters. This gives us the true measure of the

widest points in the teeth. These contact points are automatically identified and the

measurements are also automatic. The user can also change the location of the points to

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

make additional measurements of the teeth. For example, for tooth 41, the distance is 5.5

mm.

10

15

20

Figure 47B shows the teeth of Figure 47A rotated and displayed in a new orientation.

By using the camera navigation icons 766, the user can zoom in or rotate the teeth to have a

new viewpoint as desired. This enables the user to see more readily, in three dimensions, how

the teeth are oriented.

Figure 47C shows the user activating an icon 770 which uses a feature extraction

algorithm to highlight for the user the marginal ridges and indicate the distance along the

marginal ridges. The marginal ridges and associated lines connecting them across the teeth

are shown in Figure 47 C as lines 772. Other feature extraction tools are possible, including

icons for showing cusp tips and cusp fossa.

Usage of the contact points feature of Figure 47A and 47 B is shown in Figures 47 C

and 47D. In Figure 47D, the user has displayed the malocclusion teeth model 116 along

with the lines 772 indicating the marginal ridges, together with the occlusal plane 780. Note

that the lines 772 have a pronounced step relationship. In Figure 47E, the user has

interactively moved the teeth in the model 116 so that the teeth just touch the occlusal plane

780. Note that the step relationship in the lines 772 has essentially disappeared. The user can

either move a tooth individually or select a group of teeth and move them as a group. Once

the molars have been moved as shown in Figure 47D the user will typically proceed to setting

up the incisor position. To do this, the user may wish to invoke the clipping plane feature

discussed elsewhere in this document to cut through the teeth and define the upper and lower

incisor relationship in two dimensions, and then in three dimensions.

After the user has completed space management for the mandible using the tools in

the previous figures, the user proceeds to maxilla space management using the tab 790.

66

Similar screen displays as shown in the previous "Mandible Space Management" figures are

provided t o the user and thy perform space management tasks for the upper arch.

After completion of maxilla space management, the user proceeds to click on a

"Space Management (Mandible and Maxilla)" icon 792, and the screen display of Figure 48

appears. In this display, the user is provided with the same set of icons across the top of the

display, allowing the user to hide or show the upper or lower arches, gingival structure,

occlusal planes, contact points, X-ray or other image data, etc. Here, the user has selected

for display both arches in the occlused condition and the tooth model 116 is displayed. This

display includes a field 800 whereby the user can shift the midline of both the upper and

lower arches right or left by pressing the R and L buttons with the mouse. The midline then

shifts, with the amount of the shift indicated in numerical value. The display includes an

overbite icon 804, which indicates the amount of over bite in mm. The user can user the icon

to change interactively the amount of the overbite. Similarly, using icon 802 the user can

change the amount of the overjet. The tools 806 and 808 provide the length discrepancy for

both arches. Figures 48A and 48B shows similar functionality being available during

activation of the "Digital Treatment Planning Icon 900 and its associated displays; in these

figures the Target Parameter tab 902 is selected and the user is changing values for overbite

and overjet via the numerical value entry fields 910 and the results are immediately visible on

the display.

10

15

20

In Figure 49 and 49A, the user has activated the icons across the top of the display to

simultaneously display both a 2D panoramic X-ray 810 of the teeth and jaw as well as the 3D

model of the teeth, but with the teeth 116 spread out and represented in two dimensions in

registry with the 2D panorama X-ray. The teeth models that are shown in Figure 49 represent

the tooth positions in a proposed treatment. In Figure 49B, the user has unchecked the X-ray

67

icon and only the teeth are displayed. In Figures 49, 49A and 49B, the user activates the

icons to superimpose the 3D model in registry with the 2D x-ray. The user thus is able to

figure out the axis inclinations of the teeth and can revisit the selection of a reference tooth

(see Figure 27) if appropriate. Any changes in tooth position and its effect on root position

can be evaluated in the 2D/ 3D model of the panorex radiograph, CT scan, etc. plus

superimposed crowns, as in Figures 49A and 49B. From these views, the user can quickly

arrive at a proposed setup by aligning (interactively moving) the 3D teeth up and down to

find the ideal position, selecting or designing an arch form and midline, and then wrapping

the teeth around the desired arch form in three dimensions. When viewing the resulting

proposed set up, e.g., in Figure 48 and 48A, the user is again given arch length inadequacy

from a three dimensional model and proceed to modify the set up to eliminate the arch

length inadequacy using the techniques described previously. As proposed modifications are

made, the user is instantly provided with new values for arch length inadequacy, and can

immediately determine whether the proposed modifications are sufficient to remove

remaining arch length inadequacy.

10

15

20

In Figure 50 the teeth 116 that are shown are the malocclusion or original

arrangement of the teeth. Thus, Figures 49 and 50 show that the user can toggle back and

forth between initial tooth configuration and proposed treatments for the patient. Any

changes in any one-environment of module changes the values in the other environments.

In Figure 51 and Figure 51A, the user has used the navigation icons to rotate the view

shown in Figure 50. The user is provided with excellent visual aids to view how the teeth

line up with the tooth roots.

Figure 52 illustrates that the user can unclick the icon causing the X-ray to appear on

the screen and simply perform space management for both the upper and lower arches using

68

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

the virtual teeth models. As shown in Figure 52, the teeth 116 are arranged in a line and seen

from an anterior view. Figure 53 shows the view of the teeth 113 from the opposite direction.

After the user has completed the task of managing space between the virtual teeth in

the proposed arrangement; the user is able to cycle back and repeat any of the previous steps

by activating the icons on the lower left portion of the display and entering into the

appropriate displays and making further adjustments in the proposed arrangement.

The user can then access the rest of the treatment planning software, such as the

software indicated by tab 458 (Figure 21) and proceed with additional treatment planning

procedures. Finally, when they are finished, the user selects a finalized proposed treatment

plan for treating the patient. This updates the patient's prescription and is stored in the

workstation. The display of Figure 48 shows how the patient's teeth will look at the end of

treatment. Adjustments in inter-arch relationships, as shown in Figure 28 and 29, either as a

result change in overjet or overbite, or change in relationship of the jaws, are either tooth

driven or jaw driven. These actions change arch length inadequacy values. These can all be

simulated on the user interface.

10

15

20

Figure 54 shows the user entering the patient information tab 450, where the patient

has access to dental clinical examinations and dental radiographic examination. The upper

left hand portion 848 includes a slide bar 850 that allows the user to access various fields,

including demographics, patient history, examination (notes), photographs, and X-rays.

Additional fields are made available by activating the "Next" icon. The user has moved the

slide bar to X-rays. In the display 854, the user is provided with the X-rays that are currently

stored for the patient, which include a later X-ray of the face ("Lateral Ceph"). The portion

856 of the display shows the X-ray and icons for rotation of the X-ray, changing the

brightness and contrast, and displaying a mirror image of the X-ray. The user can point and

69

click in any region of interest to access dental history or access photo image databases,

radiographic images, and so forth. The navigation icons allow the user to rotate, pan, zoom

all the X-rays to see them appropriately to check for pathology, etc. Also, the user can mark

up the X-rays for making measurements in two dimensions, measuring angles, and entering

that information into the patient database.

5

10

15

20

In Figure 55, the user has navigated to a soft tissue analysis screen from the patient

information window. Here, the user is allowed to enter into the workstation via field 860

specific, numerical values for initial measurements of soft tissue, desired values for these soft

tissue parameters. In the field 862, the user is able to measure or calculate angles, e.g.,

between various teeth and various planes, angles of canting of planes that the user has

specified, etc. Again, these angles are calculated by means of display on the workstation of

patient's X-rays or other image data (including possibly scan data).

The workstation includes a database of "normal" or normative measurements for

patients of different ages, races, and sexes, for both soft tissue measurements as well as all of

the angles shown in Figure 55. The comparison thus leads the practitioner to identify

deviations from normative measurements. The display shows the normal values in the right

hand column 864 of the displays 860 and 862. Thus, as the user is designing treatment and

entering proposed or "desired" values for any of these biological parameter, the screen

display simultaneously displays the "normal" values for a patient having the same (or

approximately the same) characteristics as that of the patient.

Additional feature extraction algorithms that the workstation preferably provides

besides the marginal ridge and contact points features described previously, include

algorithms for identifying tooth cusps and fossa of the teeth. Such measurement tools are

useful in automatically performing the Bolton tooth discrepancy level and Angel

70

classification methods.

One of the unique features of the software is that the measurement features described

herein allow the practitioner to determine the Bolton tooth size discrepancy.

Bolton analysis

5

10

15

20

A method developed by W. Bolton (1958) for the evaluation of mesiodistal tooth size

discrepancies between sets of corresponding maxillary and mandibular teeth. The analysis

distinguishes between the "overall ratio," which involves all permanent teeth except the

second and third molars, and the "anterior ratio," which encompasses only the six anterior

teeth of each jaw. For this analysis it is assumed that the relatively smaller tooth material is

the correct one. A table of standard values lists the tooth width value in the opposing arch that

is ideally related to this given correct value. The difference between the ideal and actual

dental width in the arch with the excess value gives an estimate in millimeters of the severity

of tooth size discrepancy between the arches.

Tooth size discrepancy (Bolton discrepancy)

Incongruity between the sums of the mesiodistal tooth sizes of sets of corresponding

maxillary and mandibular teeth, is determined by the Bolton analysis. A discrepancy could

involve the "overall ratio" (which encompasses all permanent teeth except the second and

third molars) or the "anterior ratio" (which includes the six anterior teeth of each jaw) and is

identified as a maxillary or mandibular excess or deficiency. Only deviations that are larger

than two standard deviations are considered to be of potential clinical significance.

A tooth size discrepancy may cause difficulties in achieving an ideal overjet and

overbite or arriving at a good intercuspation during the final stages of orthodontic treatment.

Different ways to address such a problem include extraction of teeth in the arch with the

excess tooth material (usually one mandibular incisor), interproximal stripping,

71

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

compromising the angulation of some teeth so they can occupy a larger or a smaller space in

the arch, or increasing the mesiodistal tooth size in the arch with the deficiency in tooth

material (build-ups).

The present software provides measuring tools for measuring these parameters and

conducting this analysis (using the contact points algorithm described and illustrated

previously). Moreover, the workstation includes a database of normative or normal ratios for

patients. The user compares the ratio for the patient, obtained directly using the measuring

tools, and compares the result with the normative values from the database in the workstation.

The difference is displayed for the user. The result is the Bolton tooth size discrepancy and

is useful in treatment planning and allows the user to measure the total form or shape of the

teeth.

10

15

20

Another feature provided herein is the so-called "Angle classification", which is a

measure of how closely the upper and lower arches fit in an occlused condition. The

classification system is as follows.

Class I malocclusion (Neutroclusion)

A malocclusion in which the buccal groove of the mandibular first permanent molar occludes

with the mesiobuccal cusp of the maxillary first permanent molar. The term "Class I" is

sometimes used incorrectly as a synonym for normal occlusion, although in reality, it only

signifies a normal relationship of maxillary and mandibular first molars in the sagittal plane.

• Class II malocclusion (Distoclusion, Postnormal occlusion)

A malocclusion in which the buccal groove of the mandibular first permanent molar occludes

posterior (distal) to the mesiobuccal cusp of the maxillary first permanent molar. The

severity of the deviation from the Class I molar relationship usually is indicated in fractions

72

(or multiples) of the mesiodistal width of a premolar crown ("cusp" or "unit")

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

• Class II malocclusion, Division 1

A Class II malocclusion with proclined maxillary incisors, resulting in an increased

overjet

5

15

20

• Class III malocclusion (Mesiocclusion, Prenormal occlusion)

A malocclusion in which the buccal groove of the mandibular first permanent molar

occludes anterior (mesial) to the mesiobuccal cusp of the maxillary first permanent molar.

The same conventions as described above are used to indicate the severity of deviation from

a Class I molar relationship.

10 Angle classification

• "Subdivisions" (left or right) are used in asymmetric situations to indicate the side that

deviates from a Class I molar relationship.

The workstation software features measurement tools to directly make these

measurements (by measuring the distance between cusps and fossa of opposing teeth). The

results can be quantified and displayed to a user, and compared to normative values in a

database. Additionally, the values can be classified in accordance with the Angle

classification system, e.g., Class I, Class II or Class III. The resulting display of

classification is useful for interdigitation or changing the spacing between the opposing teeth.

Another feature of the software is that it allows the teeth in either or both arches to be

displayed as semi-transparent objects, which allows the user to view through the teeth to see

opposing teeth or adjacent teeth. Several possible method of providing semi-transparent teeth

is to show fewer of the points in a point cloud of teeth or fewer triangles in a mesh or triangle

73

surface representation of the teeth.

Figure 56 illustrate the patient has navigated to a cephalometric marking screen in the

patient information tab, where the user has chosen for display a lateral ceph X-ray of the

head. The user has also retracted two dimensional template teeth 870 from a library of

template teeth and superimposed the template teeth over the X-ray. The user has also

activated the icon s 872 and 873 which causes an occlusal plane 874 to appear on the display.

By activating the icons 876 in the left hand side of the display, the use can navigate to

calculations screens and associated tools, which provide the user with the ability to calculate

various parameters, a soft tissue analysis screen, landmark status, display options and other

tools. The user can move the teeth in two dimensions, both the upper and lower teeth. The

user marks the occlusal plane, and that the user is able to move the teeth relative to the

occlusal plane.

10

15

20

In Figure 57, the user has navigated to a hard tissue analysis screen, wherein the user

is provided with tools to mark various hard tissue anatomical locations on the X-ray image.

Here, the user has activated icons to compare initial and desired tooth positions using

template teeth to simulate a proposed treatment of the patient. The dark teeth 880 represent a

proposed tooth position whereas the light teeth 882 represent initial tooth positions. The user

can change the proposed tooth position by moving the dark teeth using the mouse by clicking

and dragging.

In Figure 58, the user has navigated to a screen display showing both the hard and soft

tissue, with the display 884 proving the user the tools to mark specific soft tissue locations on

the virtual patient model, including Glabella, soft tissue Naison, subnasele, mentolabial

sucus, etc. By closing out of window 884, the user accesses the window 886 where the user

is able to enter hard tissue points in a similar fashion on the display. The user is again able

to make measurements between any points that are marked on the screen and measure

74

corresponding angles.

10

15

20

The treatment planning described in figures 22-26 and 54-58 is essentially done in

two dimensions. After these initial steps are taken, the software allows the user to undergo

more specific treatment planning operations in three dimensions using the virtual 3D model

of the teeth, as described in conjunction with Figures 27-53. All changes to dental position

or bone changes can be translated into changes in soft tissue appearance using morphing

algorithms. The icons also allow for standardized views: side, planar, etc.

After the user has completed the task of managing space between the virtual teeth in

the proposed arrangement, designing the desired arch form, and arriving at a proposed tooth

arrangement or treatment plan, the user is able to cycle back and repeat any of the previous

steps by activating the icons on the lower left portion of the display and entering into the

appropriate displays and making further adjustments in the proposed arrangement.

The user can then access the rest of the treatment planning software, such as the

software indicated by tab 458 (Figure 21) and proceed with additional treatment planning

procedures. Finally, when they are finished, the user selects or saves the treatment plan. The

process can be repeated as often as desired and the screen displays are structured so that the

user can navigate anywhere in the displays at any time, and therefore repeat, as necessary, the

aligning steps, the design of the arch, enter additional patient information, access the

appliance design features and change the appliance design, etc. Moreover, as the design of

tooth finish position dictates or drives the prescription of the appliance, the present treatment

planning techniques lead directly to appliance design parameters (bracket and wire position,

75

or other treatment design, such as staged shell configuration) for treatment of the patient.

It will be appreciated that the comprehensive functionality provided by software

described herein is fully applicable to a full range of craniofacial disorders, and while the

preferred embodiment is in an orthodontic context, the invention is certainly not so limited.

It will further be noted that there may be some interdependencies between the

constraints, in other words, if the user changes one constraint, e.g., occlusal plane, other

constraints may also be affected (in 2 and 3 dimensions). Examples of such constraints

include A/P positions of the incisors and molars, intermolar width, intercanine width, amount

of overjet, amount of overbite, sagittal relation of the teeth, and lip protrusion.

Evaluation

5

10

15

20

Most of the user interface tools described herein can be used to evaluate the dentition and

the associated craniofacial structures. The dentition may be viewed in a number of states

either singularly or comparatively. There are four possible states of interest for conducting

evaluations: set-up state, initial state, response state, and final state. The set-up state is

defined as a proposed or simulated three-dimensional position of a single tooth or the

entire dentition and it reflects a state of desired or targeted occlusion. The set-up state is

arrived at by using the treatment planning software tools and features described herein (or

equivalent features). The initial state represents the non-treated malocclusion. The response

state is the current representation of the patient's orthodontic state during treatment. The

final state is defined as the state of the dentition and associated craniofacial structures at the

end of treatment. Furthermore, orthodontic appliances can be evaluated with respect to these

different states. Additionally, each of these states may be compared against normative

76

databases applicable to each state. Normative databases also exist for the appliances.

Evaluation of each state may have different goals. The objective of the treatment set-

up evaluation is to verify if the simulation meets the desired objectives of the treatment. This

is accomplished by evaluating a number of criteria: (a) boundary conditions such as arch

form, midline, reference objects such as fixed teeth (teeth to remain fixed during treatment),

etc., (b) inter arch relationships in the sagittal, vertical and transverse direction such as

overjet, overbite, etc., (c) intra-arch relationships at the individual tooth level or tooth-to-

tooth level at the first, second, or third order level. Additionally, the soft tissue structures

such as the gingival tissue or skeletal structures may be compared in association with the

setup. The patient's bite registration may be compared as well. The set-up state may be

compared to any of the previously mentioned states.

Figure 59 shows various types of treatment set-up evaluations covered under the

preferred embodiments of the present invention. Treatment set-up evaluation may be

undertaken to verify (a) treatment set-up against patient's malocclusion, (b) treatment set-up

against patient's response state, (c) appliance prescription for the treatment; while bringing to

bear in the evaluation process patient's data sets including virtual 3D models derived from in-

vivo scanning or scanning of physical model of the patient's dentition, 2D photographs, 2D

and 3D x-rays, and various combinations thereof, and practitioner notes; and treatment

related normative databases and appliance libraries including bracket libraries, arch wire

template libraries, and libraries of other types of appliances including removable shell

appliances. These treatment set-up evaluations and appliance evaluations coupled with data

sets and libraries can be done collectively in a systematic way, or individually, or in any

desired combination. Note in the present text that the evaluation of the set up state is

described comprehensively. However, the tools described herein may be effectively

employed to evaluate the additional set-up states described previously, including set-up

states representing the current situation during the course of treatment and a final position of

77

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

10

15

the teeth.

10

15

20

In preferred embodiments of the present invention, the evaluation is a combination of

user exercise of clinical judgment and software tools that are provided on the workstation.

These software tools are coded in software as instructions providing a series of predetermined

steps for guiding a user to interactively evaluate the proposed set-up. The predetermined

steps comprise steps for 1) evaluation of the proposed set-up against boundary conditions for

treatment of the patient, the boundary conditions including a least a midline and an occlusal

plane, and 2) evaluation of whether the tooth positions in both arches, and the inter-arch

relationship, of the proposed set-up are essentially ideal for treatment of the patient. As used

herein, the phrase "essentially ideal for treatment of the patient" means that the proposed set

up satisfies the practitioner's objectives for treatment of the patient, in other words, the

proposed set up meets the objectives for treatment of the patient, such as interarch

relationship, tooth position, etc., midline placement, overbite and overjet, whatever they may

be for a given patient. This invention recognizes that that there may be more than one

possible ideal setup for a given patient, as different practitioners may arrive at unique

solutions for treatment for the patient, and that for any given practitioner there may be more

than one proposed setup that meets the objectives for treatment of the patient, any one of

which would satisfy the term "essentially ideal for treatment of the patient."

The set-up evaluation may be performed by the practitioner that prepared the set up,

or it may be performed by a third party, who may be located remotely. In either situation, a

workstation is provided which includes the virtual 3D model of the proposed set-up and the

78

evaluation tools as described at length herein.

The purpose of the evaluation is several-fold:

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

a) Determine how close the set up is to the objectives identified for treatment of the

patient;

5

10

15

20

b) Determine how close the set-up is to established clinical standards of care;

c) Provide a means by which the proposed set-up can be evaluated in a multi-

practitioner environment, for example in a situation where all practitioners in the

group have to agree on the set-up. Obviously, the proposed set up is

transportable over a communications network (e.g., Internet) in the form of a file

or else accessible on a central server where others can access it using known

computer networking technology.

d) The evaluation can also serve as a guide by evaluating the course of treatment

and the eventual outcome of treatment and providing a means to measure the

difference between the actual outcome and the expected outcome. To realize

this aspect, the practitioner would need to periodically obtain updated scans of

the patient during the course of treatment and compare the current (or final)

tooth position with the expected position and use measuring tools or other

graphical devices (e.g., shading on tooth models) to quantify the amount of

variance between the actual and expected position. Tooth position data during

the course of treatment can be obtained by using the in-vivo scanner described

in the published PCT application of OraMetrix, cited previously.

e) Provide a check on the choice and design of therapeutic devices. Unless the

setup is evaluated correctly, any therapeutic device design (e.g., bracket

placement location or archwire shape) may be wrong. The evaluation provides

for dynamic change of the setup, and resulting appliance design, before the

initiation of treatment.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001 The treatment planning system described herein provides an ideal basis for performing these

evaluations. More particularly, it provides a complete and integrated workstation

environment and all the necessary processes and software tools to perform the evaluation.

Image data and actual measurements quantifying the relationship between the patient's teeth

and associated bone and facial anatomy are literally at the practitioner's fingertips. The

workstation provides ways to systematically look at the set-up, including providing

standardized views (plan, side, anterior, posterior, etc.). The workstation provides both

measurement tools and a full suite of visualization tools to both measure the proposed setup

as well as model interactively proposed changes in the set up. Complete user interaction

with the virtual patient and the proposed set up is possible, both in a single site (on one

workstation) and also in a multiple site situation where the model is shared with different

practitioners (e.g., surgeons, orthodontists, prosthodontists, etc. all treating the same patient).

The set -up evaluation function can be invoked by the user clicking on an icon or

node in a case management tree. See Figure 60, icon 1000. Highlighting this icon opens a

list of predetermined set-up evaluation steps, show in the window 1002 below the tree. The

use of predetermined set-up evaluation steps, such as the steps listed in the window 1002 (an

any substeps within any of the modules shown in window 1002), provides structure in the

evaluation review and assists the user in performing a complete, systematic and

comprehensive review of the set-up. The set-up evaluator then proceeds to execute each of

the evaluation steps that are indicated in the window 1002. While passing through the

evaluation steps, the user may evaluate whether the proposed set-up meets the boundary

conditions and the required treatment plan. To assist in this process, the user is provided with

an icon 1010, which, when activated, displays the user's notes and allows the user to review

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

80

10

15

their notes containing the treatment prescription that were previously entered when the

treatment plan was created.

If, during the evaluation, the user desires to make any modifications to the set-up, the

user can activate a "modify" icon 1004. This icon is present at every step in the set-up

evaluation process. The user is able to modify the treatment plan in any manner desired. The

changes carry over to further steps in the evaluation process. At the end of the process, the

user can choose to accept or reject the modifications, using the icons 1006 and 1008.

In preferred embodiments, a proposed set-up evaluation consists of a sequence of

instructions or process flow, realized by activation of the steps set forth in window 1002 in

Figure 60, which causes a sequence of screen displays to appear, each with tool bars and tabs

and shown in the Figures 60-90. In the preferred embodiment, two types of evaluations are

performed: firstly, were the boundary conditions met? In other words, conduct a check that

the boundary conditions in the set-up have been delivered in accordance with the proposed

prescription. This evaluation is shown in Figures 60-75 and 78 and described in more detail

below. The boundary conditions consist of two types: the boundary conditions of the

craniofascial complex as a whole, and boundary conditions of the dental complex, including

arch form, midline, occlusal plane(s) and position of one or more reference teeth. Secondly,

a methodical evaluation of the proposed set up is performed to see if the treatment is ideal for

the patient. This evaluation is performed in several discrete steps, as described in detail

20 below.

10

15

The first step, of analyzing the boundary conditions, involves evaluation of the

midline (viewing frontal photographs and panorex images), the tooth models and checking

the upper and lower midlines. Next, the occlusal planes are checked. Then the axes of the

reference tooth or teeth are checked. Then a check of any fixed teeth is made - confirm A/P

81

position of the teeth is fixed, the vertical position is held fixed, and the torques are held

fixed. The second step, of checking that the set-up is ideal, includes a check of the aesthetic

appearance of the teeth in the proposed arrangement, in various perspectives, and their

functional and occlusal relationship to the other teeth. The set up evaluation includes

evaluation or confirmation of the following: the interarch relationships: occlusion class for

both right and left sides, overjet, overbite, as well as tooth-to-tooth relationships: front intra-

arch incisor alignment for uppers and lowers (axial alignment, embrasures, contact points),

canine tooth intra-arch tooth alignment, lateral intra-arch buccal segment tooth alignment,

and tooth positions and rotations (in/out relationships, first order rotation alignment, contact

points for both upper and lower arches). The evaluation proceeds to checking the marginal

ridges, third order torque for the buccal segments and front segments (including left/right

torque symmetry).

10

15

20

The user then optionally performs a final check of the setup by comparing the

frontal photographs to the malocclusion and the set-up. After concluding the evaluation, the

user indicates their decision by activating one of the icons 1004, 1006 and 1008 shown in

Figure 60: Accept, Modify or Reject. If the user seeks to modify the set-up, they click on

Modify and then go back and use any of the tools described herein to make appropriate

changes, and then go through the evaluation checklist again. When they are finally finished,

they click on "ACCEPT" (icon 1008) and the proposed setup is saved as such in memory for

the workstation along with the date the set-up was accepted.

The tools used in the set up, appliance or other evaluation can vary depending on the

sophistication of the treatment planning software and the types of functions and features that

are provided to the user. All tools previously described for treatment planning are also

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

available for set-up, applicance, or other evaluation. In a preferred embodiment, these features include, among other things:

- The ability to superimpose X-ray or other images on the proposed set up;
- The ability to compare the proposed set-up against a variety of image data, including photographs;
- The ability to view and select two dimensional and three dimensional images in a standard view or in a view driven by the user;
- The ability to view boundary condition as shown in Figure 61 by an arch form,
- The ability to see multiple states along with boundary conditions such as shown for example in Figure 62 where malocclusion is shown by darker shaded teeth, final state is shown by light shaded teeth, and arch form is shown as a boundary condition;
- The ability to select and see a single tooth as shown in Figure 63A or a segment of teeth as shown in Figure 63B;
- The ability to view and select sectional views;
- The ability to view and select reference lines or planes;
- The ability to view axial inclinations as shown in Figure 64;
- The ability to hide certain views or elements as shown in Figures 65A 65C,
 where Figure 65A shows teeth, brackets, and arch wire; Figure 65B shows brackets
 and arch wire and hides teeth; and Figure 65C shows arch wire and hides teeth and
 brackets;
- The ability to view panoramic view;
- The ability to convert a three-dimensional image into a two dimensional image and back to three-dimensional image. For example, a patient's three-dimensional

5

10

15

dentition model can be viewed in two-dimension where subsequently teeth are

repositioned and then resulting two-dimensional dentition model is wrapped back

into the desired arch form thus recreating the 3D model which may be modified

from its original form;

5

10

15

20

- The ability to extract features of interest to the user;

- The ability to recall the practitioner's notes from computer memory and display

them along side the proposed set-up;

- The ability to superimpose the proposed set-up over the original malocclusion

model and display, using shading or color features, the changes in tooth position

between initial and final positions;

- The ability to superimpose the proposed set-up on a "graph paper" or grid feature

in which the lines on the graph as drawn in known increments (e.g., 1 mm), so as

to enable the user to quickly make estimates regarding distances between portions

of the virtual tooth models, or to use the graph paper grid for point-to-point

measurements, including displacements in 2D or 3D;

- The ability to have simultaneous multiple views on the same screen by way of

split-screen arrangement;

- The ability to be able to see transparently through one or more objects for viewing

characteristics on the otherwise hidden surface, e.g., see back side of teeth from the

front view of teeth as shown in Figure 66;

- The ability to be able to use multiple tools concurrently;

- Recall values from a normative database of orthodontic values for a large

population of orthodontic patients;

- Provided standardized views of the set-up, e.g., front, side, top, bottom, etc.

84

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

- Model effects of tooth movement to the proposed treatment condition on soft

tissue, such as gingival

The ability to manipulate at intra-arch level, inter-arch level, or on a tooth-to-tooth

basis;

5

10

15

20

- The ability to perform variety of comparisons between states or conditions; etc.

For set-up evaluation, records in the form of digital data pertaining to some or all of

the following, which are stored in the computer memory of the workstation, are available to

the practitioner or the user: the patient's clinical history, medical history, dental history, and

orthodontic history as well as 2D photographs, 2D radio graphic images, CT scans, 2D and

3D scanned images, 2D and 3D x-rays, ultrasonic scanned images, and in general, non-

invasive and sometimes invasive images, plus video, audio, and a variety of communication

records, such notes, records of office visits, patient letters or communications, etc. All

records and images are digitized. These patient data are available in the pre-processed form

as well as post-processed form. The records and images are made available through suitable

user interface icons which cause display of the images and records on the user interface. The

images are combined or superimposed to create a virtual patient 3D model that includes

surface features (soft tissue) of the patient. These digital records can be accessed by a single

practitioner or user or multiple practitioner- specialists or users through Internet. Such access

can be made simultaneously by some all or users or at different times.

Various possible evaluation steps may be taken in the set-up evaluation. One set of

steps is shown in the window 1002 in Figure 60. A slightly modified version of the series of

steps is a presently preferred embodiment and will now be described in detail. It will be

appreciated that considerable latitude may be exercised in the series of evaluation steps that

may be appropriate for a given patient, and thus the following description is offered by way

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

of example and not limitation. Moreover, the order in which the steps is not necessary

critical. Presently preferred embodiments provide the user to go through the series of steps in

any order they wish, jump from one step to another, omit steps, etc. simply by highlighting

and activating the step in the window 1002 in Figure 60.

5

10

15

20

Step 1. Check Midline, Occlusal Plane and fixed teeth

In the evaluation process, the evaluator checks for the compliance of the proposed set

up with practitioner developed boundary conditions in the craniofacial complex and in the

dental complex: the midline, the occlusal plane(s), and the position of fixed or reference teeth

(if any) whose position is not supposed to change during the course of treatment. In this

step, the display of the workstation shows a frontal 3D model of the proposed set-up and

provides the user with access (via suitable icons) to the frontal photographs of the patient's

face and head, lateral 2D ceph X-ray images, panorex images, and any other suitable 2D or

3D images. The user views these images and makes decisions about the correctness of the

location of the aesthetic occlusal planes, the location of the midline, and the placement and

position of the fixed tooth or teeth. Appropriate icon tools are provided by which the user is

able to input additional notes or comments via a dialog box, and icons by which the proposed

treatment plan can be modified or interactively studied, such as icons for treat to upper

midline, treat to lower midline, maintain midline, shift midline, shift occlusal plane, etc.

These references could be either boundary conditions or used as a starting point. The

tools used here are the radiographic examination records (X-rays), two-dimensional

photographs, and initial model of the teeth, and the various icons and user interface features

that allow the user to access these data points and interactively view them and change their

position. A key feature here is the ability to superimpose the virtual 3D teeth or 2D or 3D

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

image data showing bone and root structures of the teeth, and provide the user to freely

navigate through the models and images, view them from any orientation, and zoom in or out,

etc.

5

10

15

20

Moreover, the proposed set-up and its relationship to boundary conditions can be

observed in various combinations of hard and soft tissue, such as teeth and gingival, teeth and

lip or cheek, teeth and bone, or teeth, bone and soft tissue.

Figure 60 shows the user activating the step 1012 that leads to screen displays and

tabs 1014 that provide these functions.

Figure 67 shows the user interface in this initial boundary condition module, after

clicking on an icon that recalls from memory photographic image data for the patient.

Figure 68 shows the user interface after the user clicks on an icon that recalls X-ray

data.

Figure 69 shows the user interface viewing the virtual model of the lower arch in the

proposed set-up, with the midline plane 1016 shown. Discrepancy in tooth position from

initial to final positions are shown by usage of a shading feature, indicated in areas 1018.

In Figure 70, the user has activated icon 1020 to check the position of the midline of

the lower arch, and activated one or more of the icons 1022 to show the teeth of the lower

jaw, the discrepancy in tooth movement, slide lines, and other features.

In Figure 71, the user is evaluating the aesthetic occlusal plane 1026 of the upper arch

by activating icon 1024. The lower midline is shown by virtue of the user activating icon

1020.

In Figure 72, the user has activated an icon to hide the aesthetic occlusal plane and

instead is viewing the upper arch and the lower midline 1016, in its relationship to the upper

87

arch.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

In Figure 73, the user rotated the model to a side view by activating a side view icon

or using camera navigation icons, to thereby view the lower midline and its relation to the

upper arch in the side view.

5

10

15

20

In Figure 74, the user has activated the "graph paper" tool, wherein the upper arch is

shown superimposed on a ruling 1028 resembling graph paper, with each line in the graph

representing 1 mm. As used herein, the "graph paper" tool refers to a tool in which the

virtual model of the set-up is shown superimposed on a "graph paper" or checkerboard

background on the screen display, which the graph paper units in any appropriate scale.

In Figure 75, the user has activated icons to show the upper and lower arches to

check the functional and posterior occlusal planes as indicated in the region 1030 below the

display of the teeth.

Figure 78 shows a screen shot in which a check of the reference teeth is performed,

using the self-explanatory steps the user is prompted to perform shown in region 1030. The

model of the teeth includes axis features 1032 that show the axial orientation and inclination

of the tooth axis for all the teeth in both arches.

Step 2. Check Lower Arch Form and AP positions

At this step, the user is allowed to compare the lower arch form in the proposed set-up

with the malocclusion model. The two models can be superimposed on each other with the

change in position in the teeth indicated in a contrasting color on the user interface, as

described previously. The user is able to check the incisor positions, incisor angles, and

molar positions. The AP position and the incisal angle values in the proposed set-up are

displayed numerically on the user interface, and the user is able to compare these values with

88

the prescription that was entered previously.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

The user is also able to check the slide line, discussed previously, check the arch for

symmetry using tools such as the "graph paper" tool or with other measurement tools. The

user is also able to evaluate both visually and numerically the inter-molar width and the inter-

canine width.

5

10

15

20

The user is also able to hide/show gingival tissue by activating the gingival icons.

This allows the user to evaluate the tooth-soft tissue relationships for all teeth and at the

malocclusion contact points.

Figure 76 shows the user entering this step of the process. The lower arch is

displayed in the main window. The user is prompted to perform the steps shown in area

1040. Additional tabs 1042 are provided for modifying the arch form, modifying the inter-

arch relationship of the teeth including simulating performing extractions of teeth,

interproximal reductions, or buildup on teeth, and simulating performing mandible

displacements.

Step 3. Check upper arch form and AP positions

This step is similar to step 2, except that the evaluation is done for the upper arch

form and the teeth in the upper arch. The user compares the upper arch form in the proposed

set-up with the malocclusion model. The two models can be superimposed on each other

with the change in position in the teeth indicated in a contrasting color on the user interface,

as described previously. The user is able to check the incisor positions, incisor angles, and

molar positions. The AP position and the incisal angle values in the proposed set-up are

displayed numerically on the user interface, and the user is able to compare these values with

the prescription that was entered previously.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

The user is also able to check the slide line, discussed previously, check the arch for

symmetry using tools such as the "graph paper" tool or with other measurement tools.

The user is also able to evaluate both visually and numerically the inter-molar width

and the inter-canine width.

5

10

15

20

The user is also able to hide/show gingival tissue by activating the gingival icons.

This allows the user to evaluate the tooth-soft tissue relationships for all teeth and at the

malocclusion contact points.

Figure 77 is a screen shot of the user entering the check maxilla module, wherein the

main window displays the maxilla and the user is provided with tools and icons as indicated

for performing the operations described above.

Step 4. Check Occlusion Class, Overbite and Overjet

Next, using the right and left side views of the proposed set up, the user is directed to

check the occlusion class for the molars, premolars, and canines. The Angle classification

of the occlusion can be determined, as described above. The degree of overjet and overbite

in the proposed arrangement can be both visually observed and quantified numerically, using

the back view and lingual views of the upper and lower arches. The numerical values are

represented in a input dialog box shown below the main viewer. Clipping plane features can

be invoked which allow for viewing the cross-section of the teeth.

The analysis of the inter-arch relationship can proceed by an evaluation of teeth on a

tooth by tooth basis, by comparison of a slide line for the upper arch with the slide line of the

lower arch, evaluation of the position and axis of inclination of a reference tooth, display of

gingival tissue, or hiding gingival tissue, and evaluation of the contact points between teeth in

the same arch, the marginal ridges, cusp tips and fossa.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

Figures 79 and 80 show the screen shot of the user interface after the user has entered

the occlusion class, overjet and overbite module by clicking on the icon 1050 in window

1002. The main window shows the upper and lower arches and the user is provided with the

tools and icons as shown for performing the checks described herein. The check boxes

shown in region 1044 provide a means for prompting the user to perform specific checks and

thereby evaluate the set-up in a systematic way.

Step 5: Intra-Arch Tooth Alignment and Relation

At this step, the user is provided with suitable displays of the proposed set-up and

associated icons to check the intra-arch vertical and axial alignment of the teeth in the

maxilla and mandible. The roots of the teeth should preferably be available for viewing in

this step, hence the panorex X-ray is available and can be overlaid on the 3D tooth model. In

this step, the user may also view the proposed interproximal reduction of the teeth and space

created by such reduction. Besides the standard views of the proposed set-up virtual model,

special objects, and views are provided to help the user in checking the alignment:

For vertical alignment of the cusps and the incisal edges, superposition of the virtual

model on "graph paper" on the screen display, a display of a line connecting the cusp tips,

and panorama views.

For axial inclination, the graph paper tool, a display of numerical values of tooth axes or

lines showing the tooth axis, the occlusal or interproximal contact points, and panorama

views.

5

10

15

20

Additionally, the screen display may provide for the display of actual tooth displacement

values in three dimensions, indicating the movement of the teeth from the initial

91

malocclusion position to the proposed set-up position.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

As before, the evaluation in this step can be performed with or without display of the

gingival tissue. All standard views of the arches are made available to the user, including

clipping plane, frontal, side, panorma, etc.

5

10

15

20

Figure 74 shows the screen show of the user interface when the user has entered into this

module. The user has selected to show the teeth of the lower arch. Each tooth is shown with

the axis of the tooth indicated. The check boxes shown in region 1052 provide prompts for

assisting the user to evaluate the set-up in a systematic manner. The user activates an icon to

show the upper arch and the user proceeds to perform the same tasks for the upper arch.

Step 6. Tooth Positions, rotations and marginal ridges

In this step, the user continues to check that the set-up is ideal by evaluating the in-out

alignment of the teeth and their position, on an individual basis, the rotation of the teeth, and

the position of the marginal ridges. The user is able to activate appropriate tools to assist in

this evaluation, including graph paper, occlusal and interproximal contact points, marginal

ridge lines and embrasure line features, and possibly others. Various views of the set-up,

including frontal, side, clipping plane, panorama, etc. are made available by activating

appropriate icons on the screen display or by using camera navigation icons.

Again, as before, these evaluations can be done with a display of the gingival tissue,

or with the gingival tissue hidden.

This step can be incorporated in the menu of display screen available in the mandible

and maxilla inter-arch relationship, as shown in Figures 82-90.

Figure 82 shows the user viewing the cusp tips and marginal ridges of the lower arch.

Figure 83 shows the user viewing the cusp tips and marginal ridges of the lower arch,

but in a different view from that shown in Figure 75.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

Figure 84 shows the user viewing the cusp tips and marginal ridges of the lower arch

in a panorama layout.

Figure 85 shows the user viewing the fossa of the teeth in the lower arch, again in a

panorama view.

Figure 86 shows the user activating a measurement tool in which the width of the

crown of the teeth of the lower arch is displayed in numerical fashion above each tooth (in

mm).

5

10

15

20

Figure 87 shows a different view of the teeth of the lower arch with the measurement

tool activated.

Figure 88 shows the lower arch displayed in a different view, again with the

measurement tool activated.

Figure 89 shows the user viewing the lower arch with the contact points between the

teeth of the lower arch and the upper arch shown in color coding. In one possible

embodiment, the contacts in which a collision occurs during occlusion would be indicated in

one color, e.g., red, and areas were no collision occurs but the distance between the teeth of

the upper and lower arches is less than some threshold would be indicated with a separate

color, e.g., green. Further details on this feature are described in the co-pending patent

application of Hans Imgrund et al., Serial No. 10/137,495 filed May 1, 2002, incorporated

by reference herein.

Figure 90 shows the user entering the module for performing similar evaluations on

the maxilla. In this figure, the user has simulated extraction of a tooth.

Step 7. Third Order Torque

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606 (312) 913-0001

At this step, the torque values for each tooth are compared with the malocclusion and

checked for symmetry. The display includes a display of respective tooth displacement

values and torque values. This steps uses the graph paper took for visualizing the buccal

torque. For the front torque, the teeth can be viewed from any direction, including from the

roots. The embrasure line may be invoked.

5

10

15

20

After the user has completed the evaluation, the user enters the order information

menu by activating icon 1060, as shown in Figure 91. The user is provided with a display of

the general prescription for treatment of the patient. By scrolling down using the scroll tool

1062, the user interface displays a complete set of prescription and diagnosis information for

treatment of the patient, including appliance design data (type, model, placement location,

etc), treatment stages, etc.

Figures 92-94 show an index of possible icons on the evaluation displays of Figures

60-90. The functions provided by the icons are set forth in the index.

While the seven-step approach just described is one possible methodical approach to

setup evaluation, it will be understood that other possible approaches are possible.

Additionally, the steps need not be performed in the order recited. In particular, the user

may wish to go back and forth between the various steps and repeat them (for example,

where changes are made to the proposed set-up), jump over some steps and come back to

them later on, etc. Still other steps could be inserted into the menu of steps.

The evaluation could be carried out on the workstation in the clinic or office where the

patient is being treated. Alternatively, the proposed set-up could be transmitted over a

communications medium (such as a computer network, cable network, telephone phone line,

etc.) to a remote workstation where the evaluation occurs by a user at a remote location. The

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

(312) 913-0001

evaluation could be performed at the remote location and a revised set-up sent back to the

clinic or office where the patient is being treated for confirmation or approval. The third

party conducting the evaluation could be trained personnel at an appliance manufacturer, a

practitioner or entity providing evaluation of set-ups as a professional service, by another

physician or health care practitioner also treating the patient, or otherwise.

Still further evaluation steps could be performed, such as evaluation of the proposed

orthodontic appliance design and selection of components of the orthodontic appliance, the

proposed position of the appliance on the teeth, and stages of treatment. Additionally, the

proposed treatment may encompass hybrid appliances (combinations of brackets and wires

and other types of tooth aligning devices, such as shells).

In a related aspect, a method is provided for orthodontic evaluation for appliances of a

patient. The evaluation is facilitated by a treatment planning workstation comprising a

computing platform having a graphical user interface, a processor and a computer storage

medium. The method includes the steps of storing in the computer storage medium digitized

records including image data and virtual dentition models of patient's initial state and a

treatment set-up state, i.e., a proposed treatment position of the teeth and associated

craniofacial structures. The computer storage medium further stores a virtual appliances,

e.g., a library of virtual appliances such as brackets and archwire. The method continues with

the step of providing in the computer storage medium a set of software instructions for

graphical user interface tools for access to the digitized records for treatment set-up

evaluation of the patient, and a set of computer instructions providing a set of treatment set-

up evaluation tools for a user. The method continues with the step of using the treatment set-

up evaluation tools to evaluating appliances selected and positioned on the virtual models for

95

proposed use in achieving the treatment-set-up state.

McDonnell Boehnen Hulbert & Berghoff 300 South Wacker Drive Chicago, IL 60606

5

10

15

The appliances may include brackets and arch wires. The computer storage medium

will typically store images of the appliances (in 2D and/or 3D), include images of the

brackets and the arch wires, The appliance libraries may include templates and specifications

of the brackets and arch wires. An example screen shot is shown in Figure 95, which shows

the brackets the user has selected for treatment of the patient, including the prescription data

for the bracket, manufacturer, model, etc. When the user clicks on one of the bracket

selections, a virtual model of a bracket is displayed. The user can access a normative

database by clicking on an appropriate tab to thereby compare the prescription against data in

the normative database.

5

10

15

20

The digitized records include the patient's scan data and photographs. The method

may include the steps of evaluating the brackets against the bracket libraries, the scan data,

and the photographs.

In one embodiment the digitized records include the patient's scan data and

photographs and the method continued with the steps of evaluating the arch wires bends of 1st

order, 2nd order, and 3rd order; evaluating arch wires force systems, and evaluating the arch

wires against standard templates from the appliance libraries. The evaluation of the brackets

and arch wires may be done iteratively and/or in combination with each other. For further

details on interactive treatment planning using virtual models of orthodontic appliances, the

reader is referred to prior application serial no. 09/834,412, filed April 13, 2001.

Presently preferred and alternative embodiments of the invention have been set forth.

Variation from the preferred and alternative embodiments may be made without departure

96

from the scope and spirit of this invention.